

PHILIPS
Day-Brite
CFI

Healthcare

Multi-Med




Facilitate **comfort,** **care & control**

Philips Day-Brite / Philips CFI **Multi-Med LED patient room luminaire** balances the complex needs of patients, staff, administrators, and facility owners. With exam, ambient, and reading functions, all patient and staff needs can be met with one luminaire. The exam function provides caregivers sufficient light to keep charts and perform tasks, while the ambient and reading functions help ensure that patient needs and lighting preferences are met. Multi-Med puts comfort and control at patient and staff members' fingertips.

These human centric functions help keep patients and staff happy, which can make administrators' jobs easier too. In addition, Multi-Med's LED technology and code compliance give facility owners the efficiency and peace of mind they require. Balancing these needs is critical, and Multi-Med delivers.

Multi-Med LED patient room luminaire

- Optional integral, low voltage controller for patient control of ambient and reading functions in 2x4 models
- Acrylic and polycarbonate lens or underlay options available to help prevent incursion of particles
- Antimicrobial finish on exposed surfaces designed to inhibit the growth of microbes
- Optional LED night/chart light with integral switch to adjust light levels to 100/70/40 or 10% - night light integrated into the exam section
- Can be installed into a suspended grid, surface mounted, or in drywall ceilings



© 2018 Signify Holding. All rights reserved.
Signify reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.



Philips Lighting North America Corporation
200 Franklin Square Drive, Somerset, NJ 08873
Tel. 855-486-2216

Philips Lighting Canada Ltd.
281 Hillmount Rd, Markham, ON, Canada L6C 2S3
Tel. 800-668-9008