

A man in a pink shirt and a woman in blue scrubs are looking at a Philips Resironics CPAP machine and accessories. The man is on the left, smiling, and the woman is on the right, holding a clear plastic container with various accessories. The CPAP machine is on a table in the foreground, and a white mask is also visible. The background is a softly lit room with a lamp.

PHILIPS

RESPIRONICS

Sleep and home
respiratory care

Product and solutions **catalog**



Table of contents

- Sleep therapy 2
- Masks 28
- Oxygen 50
- Home ventilation 70
- Airway clearance 80
- RDD 84
- Sleep diagnostics 106
- Pro-Tech sensors 130
- Medical education 148
- Value-added programs 154



Sleep therapy

Natural. Comfortable. **Sleep.**

DreamStation supports long-term patient use while creating new efficiencies for your business.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.



The Dream family

At Philips Respironics, we know that good, sound sleep is vitally important at every age — to your work, your play, your physical and emotional well-being. That's why we designed the Dream family to help make treating your sleep apnea easy and comfortable. So you get the kind of sleep you need to help you feel more like yourself every day.



DreamStation

Natural. Comfortable. Sleep.

With easy-to-navigate menus, a sleek, compact design and quiet operation, DreamStation makes it simple to start, customize and continue your obstructive sleep apnea (OSA) therapy.

- Controls and interface at the front of the unit for enhanced visibility on your nightstand.
- Colorful, easy-to-use display with icon-based navigation and daily feedback.
- Most comprehensive suite of connectivity options available among the leading CPAP brands!

¹Internal assessment of 2015 competitive CPAP data comparing to ResMed Airsense10/Aircurve10 platform and Fisher & Paykel Icon series platform.



DreamWear

Like wearing nothing at all

The open-face design provides superb comfort, allowing freedom of movement and the ability to choose your most comfortable sleep position. DreamWear is designed to be unlike any mask you've experienced.

- Users rate DreamWear more comfortable, more stable, easier to use and more appealing than their prescribed mask?
- Unique in-frame air-flow and an innovative air tube connection at the top of the head.
- Light, soft and flexible material.

²2015 User Preference Questionnaire, data on file.



DreamMapper

Take charge of your sleep.

As a mobile app and online software, DreamMapper helps keep you motivated to take charge of treating your sleep apnea.

- See daily feedback on your progress.
- Set and track custom goals, and unlock achievements.
- Find guidance and informational videos to help you stay committed.
- Receive notifications about your therapy and DreamMapper compatible* equipment.

*To see which CPAP and BIPAP devices are compatible with DreamMapper, visit sleepapnea.com/mydreammapper.

Discover the whole story on the innovative, easy-to-use and comfortable Dream family.

Visit sleepapnea.com/dream

Sleepapnea.com offers tips and advice on our product, videos of patient success stories, and provides more information about the Dream family.

DreamStation

Powerful patient-driven design

DreamStation was designed with feedback from hundreds of patients to provide a simpler CPAP experience to help both new and experienced patients use therapy long-term. A suite of intuitive tools can be customized by healthcare professionals to choose the solution that's right for each patient.

DreamStation quality comes through

- 63% quieter than its closest competitor*
- Durable, easy-to-connect, top-mount tubing connections for an exceptional customer experience
- Assembled in USA

*Based on the average change in sound level during simulated patient breathing, across a range of pressure, expiratory relief, and humidifier conditions.



Acclimating patients in the early days of therapy

EZ-Start can help patients gradually adapt to their prescribed level of therapy using an innovative pressure acclimation algorithm with automatic, personalized adjustments to CPAP pressure over the first 30 days.

SmartRamp allows users to fall asleep to lower pressures than standard ramp mode. The air pressure remains low until the patient experiences an obstructive respiratory event. Then SmartRamp responds with our patented auto algorithm to resolve the patient's apneas.

Make progress every day

Each morning patients are greeted with Daily Progress Feedback — a summary of how far they've come to help motivate them to stick with their therapy. The intuitive, color display shows a simple trend of their nightly hours of use, followed by a summary of their last 30 days of good nights with more than 4 hours of use.

Sleep therapy success at your patients' fingertips


Bluetooth® is included in every DreamStation device for easy connectivity to the DreamMapper patient self-management tool. DreamMapper helps patients take an active role in their sleep management by providing daily therapy progress, helpful troubleshooting videos, and goal-setting tools. Patients can stay engaged and get the information they need fast, without having to use your staff's time to resolve common issues.

An Encore performance

DreamStation connects to our powerful EncoreAnywhere patient compliance management software. EncoreAnywhere is HIPAA compliant and makes it easy and efficient for sleep professionals to manage patient compliance and therapy.

Focus your time where it belongs

EncoreAnywhere features Adherence Profiler. It's an optional, intuitive tool that helps you identify patients according to their probability of achieving adherence, helping focus your staff's time on the patients that have the greatest likelihood of success.



Natural.
Comfortable.
Sleep.

DreamStation supports long-term patient use while creating new efficiencies for your business.



Creating efficiencies
that make a difference

More cost-effective connectivity options

Only DreamStation offers Bluetooth® Cellular and Wi-Fi connectivity options. By providing the most comprehensive suite of connectivity choices among leading CPAP brands*, you can choose the options that are most cost effective. DreamStation's innovative design has the capability to adapt as technologies continue to advance.


Making it easy to reach your patients

- Cellular provides data on a daily basis and one hour after the patient has turned the device off
- Wi-Fi delivers more data to clinicians more often and provides coverage even in weak cellular areas
- Bluetooth is included in every device for easy connectivity with DreamMapper and Alice NightOne



*Internal assessment of 2015 competitive CPAP data comparing to ResMed Aircurve10/Airsense10 platform and Fisher & Paykel Icon series platform.

Insight driven innovation



Your DreamStation CPAP Options

	US Part number	Canada Part number	HCPCS
DreamStation Auto CPAP			
DreamStation Auto CPAP (Device only)	DSX500S11	CAX500S12	E0601
DreamStation Auto CPAP w/HT humidifier and standard tube	DSX500H11	CAX500H12	E0601, E0562
DreamStation Auto CPAP w/HT humidifier and heated tube	DSX500T11	CAX500T12	E0601, E0562, A4604
DreamStation CPAP Pro			
DreamStation CPAP Pro (Device only)	DSX400S11	CAX400S12	E0601
DreamStation CPAP Pro w/HT humidifier and standard tube	DSX400H11	CAX400H12	E0601, E0562
DreamStation CPAP Pro w/HT humidifier and heated tube	DSX400T11	CAX400T12	E0601, E0562, A4604
DreamStation CPAP			
DreamStation CPAP (Device only)	DSX200S11	CAX200S12	E0601
DreamStation CPAP w/HT humidifier and standard tube	DSX200H11	CAX200H12	E0601, E0562
DreamStation CPAP w/HT humidifier and heated tube	DSX200T11	CAX200T12	E0601, E0562, A4604

Your DreamStation bi-level options

Therapy for your noncompliant, nontolerant CPAP patients

Philips Respironics revolutionized sleep therapy by introducing bi-level positive airway pressure technology to treat obstructive sleep apnea. Emulating natural breathing with pressure relief at inhalation and exhalation, this proven technology does more than provide highly effective therapy – it delivers a level of comfort that can mean the difference between compliance and non-compliance for patients having difficulty adapting to traditional CPAP therapy.

	US Part number	Canada Part number	HCPCS
DreamStation Auto BiPAP			
DreamStation Auto BiPAP (Device only)	DSX700S11	CAX700S12	E0470
DreamStation Auto BiPAP w/HT humidifier and standard tube	DSX700H11	CAX700H12	E0470, E0562
DreamStation Auto BiPAP w/HT humidifier and heated tube	DSX700T11	CAX700T12	E0470, E0562, A4604
DreamStation BiPAP Pro			
DreamStation BiPAP Pro (Device only)	DSX600S11	CAX600S12	E0470
DreamStation BiPAP Pro w/HT humidifier and standard tube	DSX600H11	CAX600H12	E0601, E0562
DreamStation BiPAP Pro w/HT humidifier and heated tube	DSX600T11	CAX600T12	E0470, E0562, A4604





DreamStation

accessories

	Part number	CA Part number
Humidification		
DreamStation humidifier	DSXH	EUXH
DreamStation water chamber	1122520	
RP-DS Humidifier Dry Box Assembly	1120668	
RP-DS Humidifier Flip Lid Seal	1120617	
RP-DS Humidifier Dry Box Inlet Seal	1120613	
Power supply		
DreamStation 80-watt power supply	1118499	
Power cord	1005894	
Tubes		
15mm Std Tube – DreamStation - RP	PR15	
15mm Heated Tube	HT15	
Filters		
Pollen filter, reusable (1 per pack)	1122446	
Ultra-fine filter, disposable (1 per pack)	1122447	
Ultra-fine filter, disposable (2 per pack)	1122518	
Ultra-fine filter, disposable (6 per pack)	1122519	
Miscellaneous		
Link module	1120293	
Nonin SpO2, assembly	1121694	
Reuseable finger sensor (clip)	936	
Oximetry finger sensor, adult, flex	953	
Flexiwrap, sens tape, adult, 25/pk	954A	
Reuseable finger sensor (soft), medium	1062050	
Reuseable finger sensor (soft), small	1062051	
Modems		
Cellular modem, oximetry capable	100600C	100610C
Cellular modem, Non-oximetry	100605C	
Wi-Fi accessory, oximetry capable	100700W	100710W
Wi-Fi accessory, Non-oximetry	100705W	

Your System One advanced therapy options

BiPAP autoSV Advanced

Therapy for your most complicated patients

The **BiPAP autoSV Advanced** system helps you support optimal patient management – for patients with Complex Sleep Apnea or Periodic Breathing (such as Cheyne–Stokes respiration). It also includes virtually all of the features and programs that make System One the most prescribed sleep therapy system by US sleep centers.¹

1. Needham & Company 3Q15 Sleep Center Survey, August 2015

2. Javaheri, et al, Clinical Application of ASV Devices, Chest 2014 Sep;146(3):858–68.

BiPAP autoSV Advanced

Targets: complex, central, mixed sleep apneas, and Cheyne–Stokes respiration patients

- Clinically proven servo ventilation algorithm
- Auto EPAP for a patient's changing needs
- Real time access to patient data via modems
- Advanced event detection
- Bi-Flex for comfort
- Encore data management and reporting
- New patient waveform reports
- Optional integrated oximetry available



	US Part number	Canada Part number	HCPCS
BiPAP autoSV Advanced, 30cm	DS960XS	CA961XS	E0471
BiPAP autoSV Advanced with humidifier, 30cm	DS960XHS	CA961XHS	E0471
BiPAP autoSV Advanced with Heated Tube humidifier, 30cm	DS960XTS	CA961XTS	E0471

BiPAP AVAPS

Automatically adapts to a patient's changing needs

Aligned with System One sleep therapy technology, the **BiPAP AVAPS** noninvasive ventilator provides the comfort of pressure ventilation and the consistent efficacy of assured tidal volume. Adapting automatically to changing patient needs, BiPAP AVAPS incorporates advanced technologies that help simplify patient care.

Part number	HCPCS	
DS1160S	BiPAP AVAPS	E0471
DS1160HS	BiPAP AVAPS with humidifier	E0471, E0562
DS1160TS	BiPAP AVAPS, Heated Tube humidifier	E0471, E0562, A4604
CA1161X (Canada)	BiPAP AVAPS	
CA1161XTS (Canada)	BiPAP AVAPS Heated Tube core package	



To learn more visit www.philips.us/bipapavaps

System One sleep therapy accessories

Description	Part number	HCPCS
Humidification		
System One 60 Series humidifier	DS6H, IN6H (Canada) E0562	
System One 60 Series Heated Tube humidifier	DS6T, IN6T (Canada) E0562, A4604 (Includes 80 watt power supply and Heated Tube)	
Water chamber (includes both 50 and 60 Series lids)	1063785	A7046
Heated Tube humidifier upgrade kit	6T-CONVKIT	A4604
Power supplies		
60 watt power supply (for H and HS models)	1091398	
80 watt power supply (for TS model only)	1091399	
Power cord 1005894		
Tubing		
System One Heated Tube	SYS1HT15	A4604
System One Performance Tube (15mm)	1074400	A7037
Performance Tubing (22mm)	1032907	A7037
Filters		
Pollen (foam) filter, reusable (1 per pack)	1035443	A7039
Pollen (foam) filter, reusable (2 per pack)	1029330	A7039
Ultra-fine filter, disposable (1 per pack)	1063096	A7038
Ultra-fine filter, disposable (2 per pack)	1035442	A7038
Ultra-fine filter, disposable (6 per pack)	1029331	A7038
Miscellaneous		
Device side panel	1099592	
Carrying case	1063857	
Data management accessories		
Data cards and accessories		
SD card (10 pack)	1063859	
SD card mailer (10 pack)	1065146	
SD card cover (attached to device)	1063858	
SD card reader	1047300	
Modems		
Wired flow modem	100200W	A9279
Wired flow modem monthly service	100201W	
Wireless flow modem	100200C, 100210C (Canada)	
Wireless flow modem monthly service	100201C	
Broadband oximetry modem and router	100200E	A9279
Broadband oximetry modem	100205E	
Modem oximetry accessories		
Oximetry sensor	Various part numbers	
Travel accessories		
Shielded DC cord	1097568	
DC battery adapter cable	532209	
Shielded 12VDC power cord system (includes shielded DC cord and battery adapter cable)	1099011	
Universal portable battery pack	1028869	



PAP accessories

to meet your patients' needs

Discover new ways to enhance your patients' therapy experience by offering a wide variety of convenient PAP accessories. Carrying these accessories can enable you to:

- Keep in contact with patients and build lasting relationships for future opportunities
- Help patients enjoy a better quality of life with accessories that allow them to conveniently travel with their sleep therapy system and easily and discreetly store it at home
- Create a retail revenue opportunity and attract additional patients to your business

DreamStation travel case

The DreamStation travel case* gives patients a convenient, compact, and discreet way to travel with their DreamStation sleep therapy system and supplies.

Convenient

The case includes a DreamStation 65W power supply that is lighter and smaller than the DreamStation 80W power supply. The included USB port powers the DreamStation** system and can be used to charge items such as a smartphone and other electronic devices during travel.

Compact

The lightweight, hard-fabric case offers an ideal combination of protection and portability. Its individual compartments accommodate the DreamStation** system and neatly organizes the supplies. An included mesh pouch holds the tubing securely in place.

Discreet

The case is designed to make traveling with sleep therapy equipment less cumbersome. A molded insert liner helps the patient thoroughly pack the PAP system components. The case can be easily attached to a roller suitcase and its zip top helps facilitate airport screening. On a plane, it easily fits in the overhead bin or under the seat.

Benefits for patients

- The hard-fabric case provides added peace of mind that the sleep therapy system will be protected during travel.
- The included power supply offers travel convenience.
- The USB port on the 65W power supply reduces the need for an additional charger.
- Molded insert liner designed to help the patient thoroughly pack the PAP device components.



Components

The case accommodates the Philips Respironics DreamStation sleep therapy system (blower only, no humidifier), the patient's mask and tubing, a 65W power supply w/ USB port, and a removable luggage strap.

Specifications

Dimensions: 10.42" (26.46 cm) x 3.94" (10 cm) x 13.78" (35 cm)

Ordering information

	Part number
Travel case, 65W power supply w/USB port	1120135
Travel case, 65W power supply w/USB port, 5/case	1124301

*DreamStation, mask, and tubing not included

**Without humidifier

PAP travel briefcase with detachable laptop bag

The PAP travel briefcase* gives patients an all-in-one carrying option that accommodates most PAP therapy systems and humidifiers, a laptop, and other small items. It makes it easy and efficient to travel with a sleep therapy system.

The travel briefcase features two bags that conveniently zip together into one carry-on. The therapy system and accessories fit snugly and securely into customizable compartments on one side of the case, while the laptop and its accessories fit securely into the ample-sized pockets of the other side. There's even space for books, magazines, and other small travel items.

Benefits for patients

- Its compartments protect the sleep therapy system and accessories during travel.
- Its zip top facilitates easy airport screening process.
- It fits within the FAA maximum carry-on size.
- The lightweight case offers easy portability.
- It's designed to easily attach to a roller suitcase.



Removable laptop bag

Detachable shoulder strap



Power supply

Tubing

DreamStation with humidifier

Components

The briefcase accommodates most PAP therapy devices and accessories and standard-size laptop computers.

Specifications

Dimensions: 16.14" w (40.99 cm) x 13.58" h (34.49 cm) x 7.87" d (19.99 cm)
Color: Black

Ordering information

PAP travel briefcase

Part number

1114784

*DreamStation, power supply and tubing not included



Mask and tubing bedside organizer

The bedside organizer gives patients a way to conveniently store their mask and tubing without cluttering a nightstand — both at home and during travel. It also keeps the components dust-free and ready for use.

The organizer features a weighted flap to allow it to hang over a nightstand or between a patient's mattress and box spring. It even features a zippered compartment for the mask and tubing, making the organizer a great travel companion.

Benefits for patients

- The organizer provides a clean, discreet place to store mask and tubing while not in use.
- Its design looks right at home in the bedroom.
- It can be used on a nightstand or in between mattress and box spring.

Specifications

Dimensions: 7.99" (20.3 cm) x 10.94" (27.8 cm) x 3.74" (9.5 cm)
Color: Light gray

Ordering information

Bedside organizer

Bedside organizer, 10/case

Part number

1122135

1124224



LiquiCell® nasal CPAP cushions

LiquiCell® is a thin, pliable cushion infused with fluid that can be applied virtually anywhere on the face, making it versatile. The cushion enhances comfort by acting as a barrier between the user's face and the mask. Gentle adhesive keeps the cushion securely in place.

Ordering information

LiquiCell nasal CPAP cushions, pack of 30

Part number

1109097

*LiquiCell is a registered trademark of Birchwood Laboratories, Inc.

Sleep therapy device comparison and selection guide

DreamStation	Modes	Comfort features		Encore data capture and reporting			Reimbursement
		Flex comfort	Resistance Control	Patterns of use	Advanced event detection	Flow waveform	
CPAP devices							
Auto CPAP	CPAP CPAP-Check Auto-Trial Auto CPAP	C-Flex C-Flex+ A-Flex	Standard	Standard	Standard	Standard	E0601
CPAP Pro	CPAP CPAP-Check Auto-Trial	C-Flex C-Flex+ A-Flex	Standard	Standard	Standard	Standard	E0601
CPAP	CPAP	C-Flex	N/A	Standard	N/A	N/A	E0601
REMstar SE	CPAP	Flex	N/A	Standard	N/A	N/A	E0601
Bi-level devices							
Auto BiPAP	CPAP Bi-level Auto bi-level	C-Flex Bi-Flex	Standard	Standard	Standard	Standard	RAD without back-up rate E0470
BiPAP Pro	CPAP Bi-level	C-Flex Bi-Flex	Standard	Standard	Standard	Standard	RAD without back-up rate E0470
Ventilatory devices							
BiPAP autoSV Advanced	Auto servo	Bi-Flex	Standard	Standard	Standard	Standard	RAD with back-up rate E0471
BiPAP AVAPS	CPAP, S, S/T, PC, T, AVAPS	Bi-Flex (only available in S mode)	Standard	Standard	Standard	N/A	RAD with back-up rate E0471
BiPAP S/T	CPAP, S, S/T	Bi-Flex	Standard	Standard	Standard	N/A	RAD with back-up rate E0471

Advanced event detection:

Patient events are detected and logged for review by the care team for optimal patient care.

Philips Respironics devices with advanced event detection detect, respond to (some models), and log the following advanced events:

- Apnea Hypopnea Index (AHI)
- Snore
- Leak
- Obstructive airway apnea (OA)
- Clear airway apnea (CA)
- Hypopnea (H)
- Respiratory-related arousal (RERA)
- Flow limitation (FL)
- Vibratory snore (VS)
- Large leak (LL)
- Periodic breathing (PB)

Any Philips Respironics device with advanced event detection gives you the ability to review the patient flow waveform report.

Matching patients with the right therapy and products

Therapy type	CPAP		Auto-CPAP		Bi-level without back-up rate		Auto Bi-level		Servo-ventilation		Bi-level with back-up rate	
Therapy system	DreamStation CPAP Pro with C-Flex+ DSX40011 (S, H, T) ¹ CAX40012 (S, H, T) ¹ (Canada)		DreamStation Auto CPAP with A-Flex DSX50011 (S, H, T) ¹ CAX50012 (S, H, T) ¹ (Canada)		DreamStation BiPAP Pro with Bi-Flex DSX60011 (S, H, T) ¹ CAX60012 (S, H, T) ¹ (Canada)		DreamStation Auto BiPAP with Bi-Flex DSX70011 (S, H, T) ¹ CAX70012 (S, H, T) ¹ (Canada)		SystemOne BiPAP autoSV Advanced 30 cm DS960 (XS, XHS, XTS) ¹ CA961 (XS, XHS, XTA) ¹ (Canada)		SystemOne BiPAP AVAPS DS 1160 (S, HS, TS) ¹ CA 1161X (S, HS, TS) ¹ BiPAP S/T DS1060 (S, HS, TS) Not available in Canada	
Target	OSA patients Short-term Auto patients ² Home sleep testing patients ²		OSA patients Short-term Auto patients ² Home sleep testing patients ²		OSA patients OSA patients non-responsive or non-tolerant of CPAP		OSA patients OSA patients non-responsive or non-tolerant of CPAP		OSA and respiratory insufficiency patients with central apnea, periodic breathing such as Cheyne-Stokes respiration (CSR), or complex apnea		OSA and respiratory insufficiency patients with hypoventilation caused by restrictive thoracic disorders, obesity hypoventilation syndrome, central apnea, and COPD Pediatric and adult use AVAPS, C-Flex, and Bi-Flex features are not cleared for pediatric use	
Algorithm	<ul style="list-style-type: none"> Corresponding Flex technology provides pressure relief Three patient-adjustable Flex settings Fixed CPAP pressure Trial period of breath-by-breath auto algorithm CPAP-Check algorithm² 		<ul style="list-style-type: none"> Flow-based proactive and reactive intelligence Delivers lowest mean average pressure levels Flow-based A-Flex pressure relief during inspiratory and expiratory phases Three patient-adjustable Flex settings 		<ul style="list-style-type: none"> Fixed IPAP and EPAP pressures Flow-based Bi-Flex pressure relief during inspiratory and expiratory phases Three patient-adjustable Flex settings 		<ul style="list-style-type: none"> Flow-based proactive and reactive intelligence Delivers lowest mean average pressure levels Flow-based Bi-Flex pressure relief during inspiratory and expiratory phases Three patient-adjustable Flex settings 		<ul style="list-style-type: none"> Dynamic, breath-by-breath adjustment of pressure support Peak flow targeted Intra-breath support Auto-titrating EPAP to alleviate the obstructive component of SDB Automatically calculated back-up rate in times of apnea 		<ul style="list-style-type: none"> Average volume assured pressure support (AVAPS) Maintains target Vt Automatically adapts to disease progressions and changing patient needs Modes: CPAP, S, S/T, T, PC Three patient-adjustable Bi-Flex settings 	
Settings	Fixed CPAP	4-20 cm	Min CPAP	4-20 cm	EPAP	4-25 cm	EPAP	4-IPAP	EPAPmin	4-25 cm ⁴	EPAP	4-30 cm
			Max CPAP	4-20 cm	IPAP	4-25 cm	IPAP	4-25 cm	EPAPmax	4-25 cm ⁴	IPAP	4-30 cm
	Ramp time	0-45 min	Ramp time	0-45 min	Ramp time	0-45 min	Ramp time	0-45 min	PSmin	0-26 cm	I-Time	0.5-3.0 sec, off
	Flex settings	1, 2, 3	Flex settings	1, 2, 3	Flex settings	1, 2, 3	Flex settings	1, 2, 3	PSmax	0-26 cm	Breath rate	1-30 BPM, off
	Adaptive Humidity	0-5	Adaptive Humidity	0-5	Adaptive Humidity	0-5	Adaptive Humidity	0-5	Max pressure	4-30 cm	Rise time	1-3
	Classic Humidity	0-5	Classic Humidity	0-5	Classic Humidity	0-5	Classic Humidity	0-5	Breath rate (Off, Auto, 4-30 BPM)		Flex settings	1,2,3
	Heated Tube	0-5	Heated Tube	0-5	Heated Tube	0-5	Heated Tube	0-5	Rise time	0-3	Target Vt ³	200-1500 ml
	Humidity level	1,2,3 ⁵	Humidity level	1,2,3	Humidity level	1,2,3	Humidity level	1,2,3	Ramp time	0-45 min	IPAPmin ³	EPAP-IPAPmax
									Flex settings	1,2,3	IPAPmax ³	4-30 cm
									Adaptive Humidity	0-5	Adaptive Humidity	0-5
									Classic Humidity	0-5	Classic Humidity	C1-C5
									Heated Tube	0-5	Heated Tube	1-5
									Humidity level	1,2,3	Humidity level	1,2,3

¹ S or XS = SD card capable
 HS or XHS = System One humidifier and SD card capable
 TS or XTS = System One Heated Tube humidifier and SD card capable
² CPAP Pro and Auto CPAP only (with Auto-Trial and CPAP-Check)
³ AVAPS
⁴ CA961 range is 4-20 cm
⁵ REMstar SE has a fixed humidity level with Heated Tubing operation



Masks

Our complete line of nasal, full-face, minimal contact, and pediatric masks reflect our commitment to innovation and meeting patient and business needs.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.

DreamWear

Like wearing nothing at all



Patients can sleep how they want with the freedom and comfort of DreamWear. Do your patients have trouble finding a comfortable sleeping position or feel restricted by CPAP tubing? DreamWear is an innovative mask created to give your patients the choice to sleep how they want.

The best of pillows and nasal masks

DreamWear's unique design offers many of the benefits of nasal and pillows masks to allow patients to have the best of both mask types. Allowing more freedom of movement and more comfort than their prescribed mask, DreamWear patients feel like they are not wearing a mask at all.*

Users told us that DreamWear gave them more freedom to move, had a better fit than their prescribed mask, and comes closest to making them feel like they don't have anything on their face during therapy.**

Key benefits

- Tubing attachment located at the top of the head to offer freedom of movement during the night
- Prevents red marks, discomfort or irritation in the nostrils or on the nose bridge**
- Users rated DreamWear more stable, more comfortable, and more appealing than their prescribed mask*
- Soft, flexible silicone frame and fabric wraps provide a soft touch on the patient's cheek

*2015 Philips User Preference Questionnaire, data on file.

**Mask does not directly contact the bridge of the nose or nostrils



Available in 4 sizes:



Part number		HCPCS
DreamWear fit pack		
1116700	DreamWear fit pack with headgear, medium frame, all cushion sizes	A7034, A7035
1116701	DreamWear fit pack without headgear, medium frame, all cushion sizes	A7034
DreamWear with headgear, small frame		
1116685	Small	A7034, A7035
1116686	Medium	A7034, A7035
1116687	Large	A7034, A7035
1116688	Medium wide	A7034, A7035
DreamWear with headgear, medium frame		
1116680	Small	A7034, A7035
1116681	Medium	A7034, A7035
1116682	Large	A7034, A7035
1116683	Medium wide	A7034, A7035
DreamWear with headgear, large frame		
1116690	Small	A7034, A7035
1116691	Medium	A7034, A7035
1116692	Large	A7034, A7035
1116693	Medium wide	A7034, A7035
DreamWear without headgear, small frame		
1116710	Small	A7034
1116711	Medium	A7034
1116712	Large	A7034
1116713	Medium wide	A7034
DreamWear without headgear, medium frame		
1116705	Small	A7034
1116706	Medium	A7034
1116707	Large	A7034
1116708	Medium wide	A7034
DreamWear without headgear, large frame		
1116715	Small	A7034
1116716	Medium	A7034
1116717	Large	A7034
1116718	Medium wide	A7034
DreamWear replacement parts		
1116740	Small nasal cushion	A7032
1116741	Medium nasal cushion	A7032
1116742	Large nasal cushion	A7032
1116743	Medium wide nasal cushion	A7032
1116750	Headgear	A7035
116745	Small frame	
1116746	Medium frame	
1116747	Large frame	
1116754	Fabric wraps	
1116752	Sizing gauge	

Amara View

High marks for performance,
no marks on the bridge of your nose

An idea whose time has come. The Amara View's incredible under-the-nose design makes it Philips' only full-face mask that won't cause red marks on the bridge of your nose. Amara View is the only Philips full-face mask that allows the user to wear glasses, as intended, to read or watch TV.*

With its amazingly innovative design, only Philips Respironics can offer the smallest and lightest choice with the widest field of vision of all leading full-face masks. Amara View covers less of a patient's face than all leading full-face masks.*

*Internal testing; data on file

Key benefits

- Innovative design prevents red marks, discomfort or irritation on the nose bridge*
- The smallest and lightest† of all leading full-face masks
- Widest field of vision of all leading full-face masks*
- Enables wearing glasses, reading, and watching TV
- Easy to use with one modular frame for all cushion sizes
- No bulky cushion and frame in front of your face and eyes



† Without quick release tube

Part number

Amara View

Part number	Description	HCPCS
1090622	Amara View small mask with headgear	A7030, A7035
1090623	Amara View medium mask with headgear	A7030, A7035
1090624	Amara View large mask with headgear	A7030, A7035
1090632	Amara View small mask without headgear	A7030
1090633	Amara View medium mask without headgear	A7030
1090634	Amara View large mask without headgear	A7030

Amara View replacement parts

1090689	Amara View quick-release tube	00000
1090692	Amara View small cushion	A7031
1090693	Amara View medium cushion	A7031
1090694	Amara View large cushion	A7031
1090696	Amara View headgear, large	A7035
1090697	Amara View headgear, standard	A7035
1090698	Amara View headgear clip, 2pk	
1090699	Amara View sizing gauge	





Pico

Pico nasal mask delivers performance for you and your patients

Practical innovation

At Philips Respironics, we created Pico with a definitive goal in mind — to create a truly appealing traditional nasal mask that makes minimal contact with the face and offers simplicity, comfort and seal. These qualities make it an ideal choice for new and existing sleep apnea patients who want minimal disruptions to their lifestyle.

Clinicians rated Pico easier to use, set up, and obtain a seal than the leading traditional nasal masks, making it ideal for your installed base optimization and new patient set-ups.*

*Data on file. Internal Marketing Survey 2013.

Key benefits

- Pico is the lightest and smallest traditional nasal mask on the market.*
- Pico is small, light, and comfortable — while providing a high level of therapy.*
- Pico comfortably fits more than 98% of patients*
- Pico offers patients a simple-to-use mask with an excellent seal that offers the freedom to watch television and read while wearing it.



Part number		HCPCS
Pico		
1104940	Pico nasal mask fitpack, small/medium, large, extra large, with headgear	A7034, A7035
1104915	Pico nasal mask, small/medium, with headgear	A7034, A7035
1104916	Pico nasal mask, large, with headgear	A7034, A7035
1104917	Pico nasal mask, extra large, with headgear	A7034, A7035
1104918	Pico nasal mask, small/medium, without headgear	A7034
1104919	Pico nasal mask, large, without headgear	A7034
1104920	Pico nasal mask, extra large, without headgear	A7034

Replacement parts		
1104934	Pico headgear	A7035
1104935	Pico headgear, reduced size	A7035
1104936	Small/medium nasal cushion	A7032
1104937	Large nasal cushion	A7032
1104938	Extra large nasal cushion	A7032
1104939	Pico sizing gauge	A7032
1104933	Pico headgear clips	
1104931	Pico quick release swivel	

Available in fitpack and 3 sizes:



Nuance

Nuance gel pillows platform

We're setting a new standard in nasal pillows masks. Philips Respironics, one of the world leaders in gel cushion technology, is bringing the comfort and sealing power of gel to nasal pillows.

As one of the latest in our groundbreaking new line of masks, Nuance gel pillows represent a choice for better sleep therapy for clinicians, home care providers and patients.

As the first and only gel pillows platform, Nuance is designed to help patients succeed by providing a more comfortable and familiar feel. With the greater comfort and seal of gel pillows, the ability to fit many different size nostrils and choice of a fabric or gel frame, Nuance is truly a great choice in nasal pillows masks.

Key benefits

- Choice between fabric or gel frame for comfort
- Non-slip headgear stays in place throughout the night
- Gel pillows provides excellent comfort and seal
- Lightweight, flexible tubing
- Causes less nostril discomfort and irritation than the leading pillows mask
- Conforms to different size nostrils



Part number		HCPCS
Nuance		
1105160	Nuance gel pillows mask, with headgear	A7034, A7035
1106194	Nuance gel pillows mask, without headgear	A7034
Nuance Pro		
1105167	Nuance Pro gel pillows mask, with headgear	A7034, A7035
1106195	Nuance Pro gel pillows mask, without headgear	A7034
Replacement parts		
1105176	Nuance headgear	A7035
1105405	Nuance RS fabric frame	
1105177	Nuance fabric frame	
1105178	Nuance Pro headgear	A7035
1105179	Nuance Pro frame	
1105173	Nuance/Nuance Pro small gel cushion	A7033
1105174	Nuance/Nuance Pro medium gel cushion	A7033
1105175	Nuance/Nuance Pro large gel cushion	A7033
1105180	Nuance/Nuance Pro swivel tube with exhalation	A7045

Nuance. A new choice. Designed around patients



Fabric frame

- All fabric frame and headgear enhances patient comfort
- Soft fabric experience for better sleep



Gel frame

- Gel-padded frame and non-slip headgear holds mask in place
- Reduces need for re-adjustment
- High-performance look

To learn more, visit <http://www.theanswer.philips.com/nuance/>



Part number		HCPCS
Wisp, fabric frame		
1094051	Wisp mask with fabric frame, with headgear	A7034, A7035
1094056	Wisp mask with fabric frame, without headgear	A7034
Wisp, clear frame		
1094050	Wisp mask with clear frame, with headgear	A7034, A7035
1094055	Wisp mask with clear frame, without headgear	A7034
Wisp youth		
1109298	Wisp youth fabric, P, S/M, L cushions with reduced size headgear	
Wisp replacement parts		
1094086	Petite cushion	A7033
1094087	Small/medium cushion	A7033
1094088	Large cushion	A7033
1112031	Extra large cushion	
1094080	Clear frame	
1094081	Fabric frame	
1094082	Wisp headgear	A7035
1094078	Wisp headgear, large size	A7035
1109307	Wisp headgear, reduced size	
1094083	Wisp headgear clips	



All three sizes in one package:



Compact design comfortably fits more than **98%** of patients



Patients said they had a **better quality of sleep** with Wisp over the leading pillows and nasal masks during a 2012 patient preference trial.



Open field of vision for a feeling of **greater independence** – patients can read and even wear their glasses

Wisp


A combination of comfort, performance and style

Developing solutions that meet the needs of a rapidly changing industry calls for fresh thinking and collaboration. At Philips Respironics we traveled around the world asking clinicians, homecare providers, and patients how we could improve their experience.

Clinicians wanted us to resolve the tradeoff between nasal and pillows masks, homecare providers wanted a solution that could reduce operational complexity by getting more patients served from one package, and patients wanted a better overall mask experience that improved their lifestyle.

Key benefits

- Compact design for a natural fit
- Soft and comfortable for a restful night's sleep
- Open field of vision for a feeling of independence
- Quick and easy to adjust for a superior seal and personalized fit
- Contemporary frame options to complement your style



TrueBlue

TrueBlue brings our best technologies together in one mask. Designed to push the limits of performance and fit, this mask sets the standard for a good night's sleep. The TrueBlue gel nasal mask with Auto Seal technology is designed to deliver a high degree of comfort, stability, and freedom of movement with minimal adjustments. This intuitive gel mask is designed for quick set-up and fitting, fewer callbacks, and easy lab titration.



Part number		HCPCS
TrueBlue with headgear		
1071800	Petite	A7034, A7035
1071801	Small	A7034, A7035
1071803	Medium	A7034, A7035
1071804	Medium wide	A7034, A7035
1071805	Large	A7034, A7035
TrueBlue without headgear		
1071806	Petite	A7034
1071807	Small	A7034
1071808	Medium	A7034
1071809	Medium wide	A7034
1071810	Large	A7034
TrueBlue with headgear DuoPack (includes two same-size cushion sets)		
1071818	Petite	A7034, A7035, A7032
1071819	Small	A7034, A7035, A7032
1071820	Medium	A7034, A7035, A7032
1071821	Medium wide	A7034, A7035, A7032
1071822	Large	A7034, A7035, A7032
TrueBlue with headgear FitPack (includes cushions in two sizes)		
1071816	Small and medium	A7034, A7035
1071817	Medium wide and large	A7034, A7035
TrueBlue replacement parts		
1071861	Petite cushion and flap	A7032
1071862	Small cushion and flap	A7032
1071863	Medium cushion and flap	A7032
1071864	Medium wide cushion and flap	A7032
1071865	Large cushion and flap	A7032
1071875	TrueBlue headgear	A7035
1071874	TrueBlue headgear, reduced size	A7035
1071888	TrueBlue gel forehead pad	
1071877	TrueBlue talon clips (pack of 2)	



Available in fitpack and 5 sizes:



ComfortGel Blue



ComfortGel Blue

Staying true to our acclaimed blue gel masks, we've retained the globally popular comfort and flexibility features you expect, while enhancing patient comfort even further. With ComfortGel Blue, we've minimized pressure points, reduced noise, and redirected exhalation air flow away from the bed partner. These and other enhancements may help to create a better patient experience and encourage greater adherence.

Part number		HCPCS
ComfortGel Blue with headgear		
1070040	Petite	A7034, A7035
1070039	Small	A7034, A7035
1070038	Medium	A7034, A7035
1070037	Large	A7034, A7035

Part number		HCPCS
ComfortGel Blue without headgear		
1070044	Petite	A7034
1070043	Small	A7034
1070042	Medium	A7034
1070041	Large	A7034

Part number		HCPCS
ComfortGel Blue with headgear DuoPack (includes two same-size cushion sets)		
1070053	Small and medium	A7034, A7035, A7032

Part number		HCPCS
ComfortGel Blue replacement parts		
1070108	Petite cushion and flap	A7032
1070107	Small cushion and flap	A7032
1070106	Medium cushion and flap	A7032
1070105	Large cushion and flap	A7032
1070104	Petite cushion	A7032
1070103	Small cushion	A7032
1070102	Medium cushion	A7032
1070101	Large cushion	A7032
1033678	Premium headgear	A7035
1059236	RS Premium headgear (smaller)	A7035
1040114	Forehead pad	
1009723	ComfortGel Blue Quick Clips (2 per pack)	

EasyLife

Part number		HCPCS
EasyLife with headgear		
1050000	Petite	A7034, A7035
1050001	Small	A7034, A7035
1050002	Medium	A7034, A7035
1050003	Medium wide	A7034, A7035
1050004	Large	A7034, A7035

Part number		HCPCS
EasyLife without headgear		
1050010	Petite	A7034
1050011	Small	A7034
1050012	Medium	A7034
1050013	Medium wide	A7034
1050014	Large	A7034

Part number		HCPCS
EasyLife replacement parts		
1050090	Petite cushion set*	A7032
1050091	Small cushion set*	A7032
1050092	Medium cushion set*	A7032
1050093	Medium wide cushion set*	A7032
1050094	Large cushion set*	A7032
1050087	EasyLife headgear	A7035
1050086	EasyLife headgear, reduced size	A7035
1050097	EasyLife Quick Clips (pack of 2)	1050097

*A cushion set includes one inner seal cushion and one outer support cushion.

Other headgear can be used with this mask including premium headgear, RS premium headgear, deluxe headgear, and small deluxe headgear.



EasyLife nasal mask

With its Auto Seal technology, the EasyLife nasal mask nearly eliminates the need for manual mask adjustments. Its innovative, lightweight design features a unique dual-cushion construction – the inner cushion creates an instant, self-adjusting seal, while the outer cushion provides comfortable support. The only manual steps necessary are intuitive headgear adjustments, making fitting as easy as possible.



Amara

Performance, appeal, and simplicity at its best

With gel and silicone options and our one “click” design, Amara delivers performance and ease of use for your patients. With the addition of a gel cushion option, Philips Respironics has made the ingeniously-simple Amara full face mask even more brilliant. Your patients can now choose the type of comfort they prefer in their cushion with gel or silicone options.

Amara offers an excellent seal with comfort and ease-of-use advantages to current traditional leading full face offerings. Both the gel and silicone options are smaller and lighter weight than the leading full face mask, and have 60% fewer parts. Amara is lightweight with a simple one “click” design and cushion options that fits over 95% of faces.

Key benefits

- Easy to fit, assemble and use
- Smallest/lightest leading traditional full-face mask
- 60% fewer parts than the leading full-face mask
- Proven to fit 95% of patients' faces

To learn more, visit <http://www.theanswer.philips.com/amara/>

Silicone

Part number		HCPCS
Amara with headgear		
1090200	Petite (reduced size frame and headgear)	A7030, A7035
1090201	Small (reduced size frame and headgear)	A7030, A7035
1090226	Medium (reduced size frame)	A7030, A7035
1090227	Large (reduced size frame)	A7030, A7035
1090267	Amara silicone fit pack (includes all four silicone cushions, reduced size frame and headgear)	
Amara without headgear		
1090210	Petite (reduced size frame)	A7030
1090211	Small (reduced size frame)	A7030
1090215	Medium (reduced size frame)	A7030
1090216	Large (reduced size frame)	A7030
Amara replacement parts		
1090290	Petite cushion	A7031
1090292	Small cushion	A7031
1090293	Medium cushion	A7031
1090294	Large cushion	A7031
1090296	Amara headgear, reduced size	A7035
1090297	Amara headgear, standard	A7035
1040114	Forehead pad	
1071877	Talon headgear clips	
1081899	Entrainment valve swivel	



Available in 4 sizes:



Gel

Part number		HCPCS
Amara gel with headgear		
1090400	Petite (reduced size frame and headgear)	A7030, A7035
1090401	Small (reduced size frame and headgear)	A7030, A7035
1090405	Medium (reduced size frame)	A7030, A7035
1090406	Large (reduced size frame)	A7030, A7035
1090467	Amara gel fit pack (includes all four gel cushions, reduced size frame and headgear)	
Amara gel without headgear		
1090410	Petite (reduced size frame)	A7030
1090411	Small (reduced size frame)	A7030
1090415	Medium (reduced size frame)	A7030
1090416	Large (reduced size frame)	A7030
Amara gel replacement parts		
1090490	Petite cushion	A7031
1090492	Small cushion	A7031
1090493	Medium cushion	A7031
1090494	Large cushion	A7031
1090296	Amara headgear, reduced size	A7035
1090297	Amara headgear, standard	A7035
1040114	Forehead pad	
1071877	Talon headgear clips	
1081899	Entrainment valve swivel	
Amara starter kits (includes gel mask with same size silicone cushion)		
1106618	Petite	A7030, A7035
1106619	Small	A7030, A7035
1106620	Medium	A7030, A7035
1106621	Large	A7030, A7035



Available in 4 sizes:



ComfortGel Blue Full



ComfortGel Blue Full full-face mask

The design combines from very familiar – and popular – components. Our comfortable forehead cushion and enhanced StabilitySelector combine with our popular blue gel replaceable cushion to create a level of comfort and stability no other full-face mask can approach. The unique sure seal technology (SST) provides a membrane-thin layer over the gel cushion to follow the contours of the face better, allowing wearers to quickly achieve the comfortable seal they need.

Other headgear can be used with this mask including premium headgear, RS premium headgear, deluxe headgear, and small deluxe headgear.

Part number		HCPCS
1040135	ComfortGel Blue Full mask with premium headgear, small	A7030, A7035
1040136	ComfortGel Blue Full mask with premium headgear, medium	A7030, A7035
1040137	ComfortGel Blue Full mask with premium headgear, large	A7030, A7035
1062004*	ComfortGel Blue Full mask with premium headgear, extra large	A7030, A7035
1002854	ComfortGel Blue Full mask DuoPack with premium headgear, small plus small cushion	A7030, A7035, A7031
1002855	ComfortGel Blue Full mask DuoPack with premium headgear, medium plus medium cushion	A7030, A7035, A7031
1002856	ComfortGel Blue Full mask DuoPack with premium headgear, large plus large cushion	A7030, A7035, A7031
1062005*	ComfortGel Blue Full mask DuoPack with premium headgear, extra large plus extra large cushion	A7030, A7035, A7031

Accessories		HCPCS
1040568	ComfortGel Blue Full mask without headgear, small	A7030
1040569	ComfortGel Blue Full mask without headgear, medium	A7030
1040570	ComfortGel Blue Full mask without headgear, large	A7030
1062003*	ComfortGel Blue Full mask without headgear, extra large	A7030
1051032	ComfortGel Blue Full mask without headgear, DuoPack, small plus small cushion	A7030, A7031
1051033	ComfortGel Blue Full mask without headgear, DuoPack, medium plus medium cushion	A7030, A7031
1051034	ComfortGel Blue Full mask without headgear, DuoPack, large plus large cushion	A7030, A7031
1040798	ComfortGel Blue Full cushion, small	A7031
1040799	ComfortGel Blue Full cushion, medium	A7031
1040800	ComfortGel Blue Full cushion, large	A7031
1062030*	ComfortGel Blue Full cushion, extra large	A7031
1040795	SST flap only, small	A9999
1040796	SST flap only, medium	A9999
1040797	SST flap only, large	A9999
1062031*	SST flap only, extra large	A9999
1044861	ComfortGel Blue Full cushion and SST flap, small	A7031
1044862	ComfortGel Blue Full cushion and SST flap, medium	A7031
1044863	ComfortGel Blue Full cushion and SST flap, large	A7031
1062032*	ComfortGel Blue Full cushion and SST flap, extra large	A7031
1040138	ComfortGel Blue Full premium headgear	A7035
1040114	Premium forehead pad	
1049398	Entrainment valve swivel	
1040139	Ball-and-socket Quick Clips, blue, 2 per pk.	
1010872	Port caps, 5 per pk.	
1049507	Sizing gauge	

FitLife

Part number		HCPCS
FitLife with headgear		
1060801	Small	A7030, A7035
1060802	Large	A7030, A7035
1089991	Extra large	A7030, A7035

FitLife without headgear		
1072530	Small	A7030
1072531	Large	A7030
1089992	Extra large	A7030

FitLife replacement parts		
1060807	FitLife headgear, small	A7035
1060808	FitLife headgear, large	A7035
1089990	FitLife headgear, extra large	A7035
1060913	FitLife snap clips, 2 per pk.	
7041	Mask swivel	
1060809	Entrainment elbow	A7045



FitLife

A solution for patients who have difficulty wearing masks or who simply prefer a clear line of sight. FitLife addresses a number of patient issues. Facial hair, dentures, nasal bridge breakdown, facial irregularities, and claustrophobia are issues that can greatly impact the ability to establish effective therapy and comfort. FitLife quickly establishes a robust seal around the less pressure-sensitive perimeter of the face, while maintaining an unrestricted view.



OptiLife nasal mask

Features a headgear design with an integrated chin support band that allows it to be put on and adjusted with little effort – and no buckles.

Part number		HCPCS
OptiLife		
1036800	OptiLife FitPack with headgear and chin support band size P, S, M, L pillows cushions	A7034, A7035
1036801	OptiLife FitPack/DuoPack with headgear and chin support band size P, S, M pillows cushions and S, M CradleCushions	A7034, A7035, A7032
1036802	OptiLife FitPack/DuoPack with headgear and chin support band size M, L pillows cushions and M, L, LN CradleCushions	A7034, A7035, A7032
1036818	OptiLife mask sub-assembly with headgear and chin support band No cushions	

Accessories		HCPCS
1036819	OptiLife sub-assembly No cushions or headgear	
1036804	OptiLife FitPack, mask only size P, S, M, L pillows cushions	A7034, A7035
1036838	OptiLife pillows cushion, petite	A7033
1036839	OptiLife pillows cushion, small	A7033
1036840	OptiLife pillows cushion, medium	A7033
1036841	OptiLife pillows cushion, large	A7033
1036846	OptiLife CradleCushion, small	A7032
1036847	OptiLife CradleCushion, medium	A7032
1036848	OptiLife CradleCushion, large narrow	A7032
1036849	OptiLife CradleCushion, large	A7032
1036850	Headgear Includes chin support band	A7035
1039967	Tube management clips, 2 per pk.	
1032007	Respironics pillows sizing gauges, 5 per pk.	
1031154	Respironics pillows sizing gauges, 20 per pk.	

*Requires addition of mask cushion for HCPCS A7034.

Accessories (cont.)

302175	Chin strap
1012911	Premium chin strap
302425	Deluxe chin strap
1032907	Lightweight CPAP tubing (6 ft.)
7212	Standard CPAP tubing (10 ft.)
1036425	Sizing gauge ring

Codes for nasal masks

A7034	Nasal mask only
A7035	Headgear
A7032	Nasal cushion replacement

Codes for full-face masks

A7030	Full-face mask only
A7035	Headgear
A7031	Full-face cushion replacement

Codes for nasal pillows masks

A7034	Nasal mask only
A7035	Headgear
A7033	Nasal pillow replacement
A7032	Nasal cushion replacement (OptiLife CradleCushion)
A7045	Exhalation port with swivel (GoLife PN 1078407)



CradleCushion



Pillows cushion

Oxygen

PHILIPS
RESPIRONICS

Philips Respironics Right Fit oxygen therapy products are intended to give patients the freedom they want while helping homecare providers reduce expenses and make the most of shrinking reimbursements.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.

SimplyGo Mini

A reliable, powerful choice in lightweight portable oxygen concentrators.

We built SimplyGo Mini to be tough. During development, SimplyGo Mini was rigorously tested in real world environments including impacts, vibrations, and temperatures to deliver reliable performance day in and day out.

Delivering more when you're on the go

From a walk in the park to an overnight trip, SimplyGo Mini may help you remain active and provide you with flexibility while on the go.

Part number

HCPCS

SimplyGo Mini System

System with standard battery - U.S.

1113601	Includes SimplyGo Mini portable oxygen concentrator, carrying case, one standard rechargeable Lithium Ion battery, DC and AC power cords, accessory bag, and user manual	E1390/E1392
---------	--	-------------

System with extended battery - U.S.

1113602	Includes SimplyGo Mini portable oxygen concentrator, carrying case, one extended rechargeable Lithium Ion battery, DC and AC power cords, accessory bag, and user manual	E1390/E1392
---------	--	-------------

SimplyGo Mini with standard battery bundles

1113601/1020000	Bundle w/EverFlo	E1390/E1392
1113601/1020001	Bundle w/EverFlo OPI	E1390/E1392
1113601/1020014	Bundle w/EverFlo Q	E1390/E1392
1113601/1020015	Bundle w/EverFlo Q OPI	E1390/E1392

SimplyGo Mini with extended battery bundles

1113602/1020000	Bundle w/EverFlo	E1390/E1392
1113602/1020001	Bundle w/EverFlo OPI	E1390/E1392
1113602/1020014	Bundle w/EverFlo Q	E1390/E1392
1113602/1020015	Bundle w/EverFlo Q OPI	E1390/E1392



The manufacturer of this POC has determined this device conforms to all applicable FAA requirements for POC carriage and use on board aircraft.



Accessories

1116816	Standard lithium ion battery
1116817	Extended lithium Ion battery
1116819	SimplyGo Mini DC power supply - U.S.
1116818	SimplyGo Mini AC power supply - U.S.
1082664	SimplyGo Mini AC power cord - U.S.
1116830	SimplyGo Mini external battery charger - U.S.
1118724	User manual

Bags - Brown

1119897	SimplyGo Mini carry bag and strap, brown
1116824	SimplyGo Mini accessory bag, brown
1116836	SimplyGo Mini backpack, brown

Bags - Black

1119928	SimplyGo Mini carry bag and strap, black
1119893	SimplyGo Mini accessory bag, black
1119894	SimplyGo Mini backpack, black

Warranties

1120121	SimplyGo Mini 3-year system, 1-year sieve warranty
1122070	SimplyGo Mini 3-year system, 3-year sieve warranty
1120123	SimplyGo Mini 2-year sieve warranty extension

SimplyGo



Part number		HCPCS
DreamWear fit pack		
1068987	SimplyGo portable oxygen concentrator	E1390, E1392
1068787, 1020000	SimplyGo portable oxygen concentrator with EverFlo non-OPI oxygen concentrator	E1390, E1392
1068787, 1020001	EverFlo OPI oxygen concentrator	E1390, E1392
1068787, 1020015	EverFlo Q non-OPI oxygen concentrator	E1390, E1392
1068787, 1020014	SimplyGo portable oxygen concentrator with EverFlo Q with OPI-oxygen concentrator	E1390, E1392
Accessories		
1082663	SimplyGo carry case	
1083692	SimplyGo DC power supply	
1082661	SimplyGo AC power supply	
1082664	SimplyGo AC power cord	
1083693	SimplyGo airline power cord	
1083696	SimplyGo accessory bag	
1106598	SimplyGo shoulder strap	
1083699	SimplyGo humidifier pouch	
1101172	SimplyGo humidifier connecting tube kit	
1101602	SimplyGo humidifier connecting tube and pouch	
1082662	SimplyGo lithium ion battery with handle	
900-103	Battery charger	
1074885	Mobile cart	
1083700	User manual	
Warranties		
1097623	SimplyGo two-year total warranty	
1097577	SimplyGo three-year total warranty	
1097778	SimplyGo five-year total warranty	

SimplyGo portable oxygen concentrator

SimplyGo is the only portable oxygen concentrator (POC) to offer 2x more oxygen than other POCs weighing 10 pounds or less. That means you can manage the therapy and lifestyle needs of nearly all of your patients with just one POC. Designed to be rugged and easy to use, SimplyGo is truly a new breed of lightweight and attractive portable concentrator.

Now one POC meets the portable oxygen needs of nearly all of your patients.

SimplyGo is Philips first oxygen concentrator that offers the benefits of a continuous flow mode in a lightweight, easy-to-carry design of 10 lbs or less.

No need to compromise between an effective oxygen therapy output for your patients and the usual weight of the device. SimplyGo offers a solution that fulfills the needs of most of your oxygen-therapy patients. With SimplyGo, patients have the flexibility to be away for the day, for work, for a weekend or a long holiday.

To learn more visit www.simplygo.respironics.com

Key benefits

- Attractive, sleek design
- Easy-to-remove battery
- Intuitive, easy-to-read screen
- Convenient carrying case
- Handy fold-up cart
- Detachable accessory bag

Simple-to-use touch-pad interface

Simple-to-select oxygen delivery mode

Lightweight Simple-to-carry anywhere, anytime!



The United States Federal Aviation Administration (FAA) has approved the SimplyGo for in-flight use on commercial airliners.





SimplyFlo

stationary oxygen concentrator

Life made simple

For nocturnal patients who may be refusing or resisting therapy because of the size, weight and lifestyle disruptions that can accompany the use of traditional concentrators, SimplyFlo weighs less than nine pounds and is designed to enhance adherence. SimplyFlo is one of the smallest and lightest continuous flow stationary oxygen concentrators available. It operates on external power only and has no internal battery. An optional travel accessory case allows oxygen patients to easily transport SimplyFlo between locations.

Key benefits

- Weighs 8.5 pounds
- Provides .5 to 2 LPM continuous flow output
- Operates on AC power (no internal battery)
- Uses less than 1/3 the electricity of traditional stationary concentrators
- Simplifies transport for travel or within the home
- Helps nocturnal oxygen users maintain an active lifestyle during the day

Part number

11076773 SimplyFlo oxygen concentrator (includes concentrator, AC power cord and user manual)

HCPCS

E1390

Accessories

1076830	SimplyFlo travel case
1082664	AC power cord
1082661	AC power supply
1076809	User manual
1076829	Humidifier bottle holder
1101172	Humidifier connecting tube kit

Warranties

1097623	SimplyFlo two-year total warranty
1097577	SimplyFlo three-year total warranty



Another solution for nocturnal patients

The optional travel accessory case allows patients to easily transport SimplyFlo between locations.

Once again, Philips Respironics is bringing innovative solutions to your business and your patients. With SimplyFlo, Philips has created the smallest and lightest stationary oxygen concentrator.

SimplyFlo represents a comprehensive option for nocturnal patients who may be refusing or resisting therapy because of the size, weight and lifestyle disruptions that can accompany the use of traditional concentrators.

Designed to enhance adherence

For many patients new to therapy who only require nocturnal treatment, the intimidating size, weight, and sound of a traditional oxygen concentrator can be a large factor in their rejection of therapy.

Until SimplyFlo, there has not been a more transportable stationary concentrator designed to help nocturnal oxygen patients.



Concentrators

EverFlo concentrator

Delivers a great performance and demands so little in return. It's small, quiet, and designed to be easier and less expensive to operate. Plus, EverFlo reduces scheduled maintenance visits. So now you can control oxygen therapy costs and deliver all the features you and your patients want with one 31-pound concentrator that's bound to be the star of your business.



Part number		HCPCS
1020000	EverFlo oxygen concentrator	E1390
1020001	EverFlo oxygen concentrator with oxygen percentage indicator (OPI)	E1390
1039362	EverFlo oxygen concentrator, transfill	E1390
1039363	EverFlo oxygen concentrator with OPI, transfill	E1390

Units include a patient manual.

Accessories

Filters

1038831	Compressor intake filter
520	Micro disk filter

Tubing

1039642	Humidifier connector tube
---------	---------------------------

Miscellaneous

H537-1	Stand-alone pediatric flow meter
H644	Integrated pediatric flow meter
528	Flow meter locking kit
H649	Caster brake kit
1038571	Warranty, two-year with compressor

EverFlo Q ultra quiet 40 decibel stationary oxygen concentrator

Our compact, 31-pound EverFlo Q is the first ultra-quiet concentrator available at an affordable price. The new EverFlo Q combines the small size, light weight, and low maintenance of our popular EverFlo oxygen concentrator, with the added value of a quiet sound level. With EverFlo Q, you get the benefits of a quiet concentrator without the burden of paying a premium which is good for you and your patients.



Part number		HCPCS
1020015	EverFlo Q oxygen concentrator	E1390
1020014	EverFlo Q oxygen concentrator with OPI	E1390
1039365	EverFlo Q oxygen concentrator, transfill	E1390
1039364	EverFlo Q oxygen concentrator with OPI, transfill	E1390

*Units include a patient manual.

Accessories

Filters

1038831	Compressor intake filter
520	Micro disk filter

Tubing

1039642	Humidifier connector tube
---------	---------------------------

Miscellaneous

H537-1	Stand-alone pediatric flow meter
H644	Integrated pediatric flow meter
528	Flow meter locking kit
H649	Caster brake kit
1038571	Warranty, two-year with compressor





UltraFill

Home oxygen filling system

UltraFill is an advanced home oxygen system that combines a stationary oxygen concentrator, filling station, and high-capacity cylinders to meet the needs of a wide range of oxygen patients, including those who are highly active or require continuous flow oxygen. UltraFill also has the needs of your business in mind with features intended to save you money by using existing inventory and lowering your operating costs.

Patient-friendly design

- One-handed cylinder connect/disconnect
- Patient can stop and restart filling at their convenience
- Low noise level: 44dBa
- Cylinders can be topped off to make partial fills
- Low power consumption: 110W
- Fill level indicator monitors cylinder fill level



Part number		HCPCS
1057100	UltraFill station 120V	K0738
1057101	UltraFill system with 2 cylinders, 2,000 psi and a cylinder bag, 120V, Canada	K0738
1074932	UltraFill system with 2 cylinders, 2,000 psi and a cylinder bag, 120V 3,000 psi cylinders available for additional charge.	K0738
1039362	EverFlo 120V, transfill	E1390
1039363	EverFlo 120V OPI, transfill	E1390
1039365	EverFlo Q 120V, transfill	E1390
1039364	EverFlo Q 120V OPI, transfill	E1390
M10600-10	Millennium 120V, transfill	E1390 (with QF or QG modifier)
M10605-10	Millennium 120V OPI, transfill	E1390 (with QF or QG modifier)

Accessories

Cylinders

1065705	M4 cylinder, US, 2000 psi
1065706	M6 cylinder, US, 2000 psi
1065707	M9 cylinder, US, 2000 psi
1065708	MD15 cylinder, US, 2000 psi
1065709	ME24 cylinder, US, 2000 psi
1065710	M4 cylinder, Canada, 2000 psi
1065711	M6 cylinder, Canada, 2000 psi
1065712	M9 cylinder, Canada, 2000 psi
1065713	MD15 cylinder, Canada, 2000 psi
1065714	ME24 cylinder, Canada, 2000 psi
1065715	MB08 cylinder, US, 3000 psi
1065716	MC13 cylinder, US, 3000 psi
1065717	MD22 cylinder, US, 3000 psi
1065718	ME36 cylinder, US, 3000 psi
1065719	MB08 cylinder, Canada, 3000 psi
1065720	MC13 cylinder, Canada, 3000 psi
1065721	MD22 cylinder, Canada, 3000 psi
1065722	ME36 cylinder, Canada, 3000 psi
1113356	MB08 3000 PSI Cylinder with OCD
1113355	MC13 3000 PSI Cylinder with OCD
1113354	MD22 3000 PSI Cylinder with OCD
1113353	ME36 3000 PSI Cylinder with OCD
1116898	MB08 Cylinder with OCD, CANADA
1116899	MC13 Cylinder with OCD, CANADA
1116900	MD22 Cylinder with OCD, CANADA
1116901	ME36 Cylinder with OCD, CANADA

Miscellaneous

1065697	Chad Therapeutics, Bonsai oxygen conserving device
1065698	Chad Therapeutics, 3100 pediatric regulator
1065699	Chad Therapeutics, 4800 regulator
1065700	Chad Therapeutics, 3100 regulator
1065701	Small cylinder bag (M9/MC13)
1065702	Medium cylinder bag (M4/M6/MB08)
1065703	Large cylinder bag (MD15/MD22)
1065724	Cylinder cart (MD15/ME24/MD22/ME36)
1065696	UltraFill transfer tubing



Millennium M10

Concentrator

Delivers up to 10 LPM of oxygen reducing the delivery costs associated with 5 LPM and above oxygen patients. The M10 is engineered to reduce the cost of providing oxygen at higher liter flow and has fewer parts than any other concentrator on the market. Fewer parts, combined with great reliability, help you to reduce expenditures on accessories and unnecessary inventory.

Part number		HCPCS
M10600	Millennium 10-liter oxygen M10 concentrator	E1390
M10600-80	Millennium 10-liter oxygen M10 concentrator, Canada	
M10605	Millennium 10-liter oxygen M10 concentrator with OPI, Canada	E1390
M10605-80	Millennium 10-liter oxygen M10 concentrator with OPI, Canada	
M10600-10	Millennium 10-liter oxygen concentrator, transfill	E1390
M10605-10	Millennium 10-liter oxygen concentrator with OPI, transfill	E1390

Both units include a patient manual.

Accessories		
1014653	High flow starter kit Includes cannula, 25' tubing, connector	A4615*, A4616*
1014654	High flow nasal cannula	A4615*
1017796	Inlet filter replacement kit	
1017768	High flow oxygen tubing, 25'	A4616*
1014655	Oxygen tubing connector	

*Accessories are not separately billable to the Medicare program, with the exception of replacement circuits. Consult your local Medicare DME MAC local coverage determination (LCD) for details. Commercial insurance plans may provide separate coverage for accessories. Consult insurance company policy manual for more information.

Oximeters



920M Plus hand-held oximeter with memory

The 920M Plus hand-held pulse oximeter provides on-the-spot oximetry readings, and is designed to be easy to use for both healthcare professionals and patients. The 920M Plus features overnight recordings (72 hours of memory), is compatible with virtually all Respironics sensors, uses disposable or rechargeable batteries, and interfaces with Profox and nVISION software.

Part number

920MP	920M Plus hand-held oximeter with memory Includes 936 adult finger clip sensor, four AA batteries, and a manual.
1028009	2500A hand-held oximeter downloadable with alarms

Accessories

Software

H920M-40	Profox software Download cable sold separately.
1019441	nVision software Download cable sold separately.

Sensors

1062051	Soft sensor, small 1M
1062050	Soft sensor, medium 1M
1062049	Soft sensor, large 1M

Miscellaneous

920MPC-NA	920MP rechargeable battery system, wall mount Includes 920MPC charger stand, NiMH rechargeable battery, power supply, and an operator manual.
1053371	Download cable
920MPC	Charger stand only, replacement
1007355	920MP NiMH rechargeable battery
1007332	Carrying case

See the replacement parts catalog for additional replacement parts.

2500A hand-held oximeter with memory and alarms

The 2500 hand-held oximeter is a small, versatile digital pulse oximeter for accurate assessment of blood oxygen saturation and pulse rate. It is ideal for ambulatory monitoring and overnight recordings. The 2500A has two-button operation, is compact, has three power options, and features visual and audible alarms. Compatible with nVision and Profox software.



Part number

920	Hand-held oximeter
920M	Hand-held oximeter with memory

Both include a 936 adult finger clip sensor, six AA batteries, and a manual.

Accessories

Hand-held oximeter

100-0433-00	Battery door
927	Extension cable, 3'
1008130	Extension cable, 10'
924	Manual
1007332	Carrying case

Hand-held oximeter with memory

H920M-40	Profox software Download cable sold separately.
1053371	Download cable
1019441	nVision software Download cable sold separately.
100-0433-00	Battery door

See the replacement parts catalog for additional replacement parts.

Part number

1062232	GO ₂ finger oximeter
---------	---------------------------------

Accessories

1062965	Carrying case
1062966	Lanyard

Part number

1078606	WristOx ₂ 3150, pulse oximeter with standard accessories (reusable soft sensor, wrist bands, batteries, CD)
---------	--

Accessories

1078747	Medium soft sensor, 8 pin (12"/0.3 m)
1078749	Adult flex sensor, 8 pin (12"/0.3 m)
1079107	Adult finger clip sensor, 8 pin (12"/0.3 m)
1084509	Sensor adapter, WristOx ₂ , 3150 (8 to 9 pin)
1078748	Sensor adapter, WristOx 3100 (9 to 8 pin)
1078750	USB PC download cable for WristOx2
1079118	Wrist band, reusable (10")
1091490	WristOx ₂ , wrist band (8")
1091489	WristOx ₂ , wristband (6")



920 and 920M hand-held oximeters

The 920 and 920M oximeters are small, easy-to-use hand-held devices that deliver on-the-spot oximetry readings; perfect for home use. The 920M has the added capability of storing up to 18 hours of trend memory. Both oximeters are intended for spot check use and can be connected to an external printer.



GO₂ pulse oximeter

The quality construction finger oximeter allows your patients to take greater control of their lives while creating a new opportunity for your business.



WristOx₂ ambulatory digital pulse oximeter

WristOx₂ is a wrist-worn pulse oximeter that has been developed to deliver readings and help reduce the difficulties that can be associated with overnight oximetry studies.

Oxygen accessories

Part number	
Cannulas	
1025655	Cannula, 4'. For use with ePOD.
1014653	High flow starter kit Includes cannula, 25' tubing, connector. For use with M10, M10600, and M10605 oxygen concentrators.
1014654	High flow nasal cannula For use with M10, M10600, and M10605 oxygen concentrators.
1017768	High flow oxygen tubing, 25' For use with M10, M10600, and M10605 oxygen concentrators.
1014655	Oxygen tubing connector For use with M10, M10600, and M10605 oxygen concentrators.
Carry Bags/Accessory Bags	
1119897	SimplyGo Mini Carry Bag and Strap, Brown
1116824	SimplyGo Mini Accessory Bag, Brown
1116836	SimplyGo Mini Backpack, Brown
1119928	SimplyGo Mini Carry Bag and Strap, Black
1119893	SimplyGo Mini Accessory Bag, Black
1119894	SimplyGo Mini Backpack, Black
1082663	SimplyGo carrying case
1083696	SimplyGO Accessory bag
1119897	SimplyGo Mini Carry Bag and Strap, Brown
1116824	SimplyGo Mini Accessory Bag, Brown
1116836	SimplyGo Mini Backpack, Brown
1119928	SimplyGo Mini Carry Bag and Strap, Black
1119893	SimplyGo Mini Accessory Bag, Black
1119894	SimplyGo Mini Backpack, Black
900-107-02	Black carrying case, EverGo oxygen concentrator
900-250	Accessory case, EverGo oxygen concentrator
1018019	GoLox carrying case
1062965	GO ₂ carrying case
1019439	WristOx ambulatory digital pulse oximeter carrying case
1007332	Oximeter carrying case, 920, 920M, and 920MP oximeters.

Part number	
Carts	
1065724	UltraFill Cylinder Cart (MD15/MD22/ME24/ME36)
1074885	SimplyGo Mobile cart
900-100	Mobile cart For use with the EverGo and SimplyGo oxygen concentrators.
Extension cables	
927	Extension cable, 3' For use with 920, 920M, and 930 oximeters.
1008130	Extension cable, 10' For use with 920 and 920M oximeters.
Filters	
1038831	Compressor intake filter For use with EverFlo oxygen concentrators.
900-150	Air inlet filters, 6 per pk. For use with EverGo oxygen concentrators.
520	Micro disk filterE For use with EverFlo oxygen concentrators.
1017796	Inlet filter replacement kit For use with M10, M10600, and M10605 oxygen concentrators.
1015581	Kit, grommet, 10 pk.
H619	Air inlet filter replacement kits, 6 per pk. For use with M600 and M605 oxygen concentrators.
H620	Inlet filter replacement kits, 6 per pk. For use with 600, 605, M600, and M605 oxygen concentrators.
H621	Microdisk filter replacement kits, 6 per pk. For use with all Millennium oxygen concentrators.
H652	Filter cover for compressor kit For use with 600, 605, M600, and M605 oxygen concentrators.
H653	Filter felt for compressor kits, 6 per pk. For use with 600, 605, M600, and M605 oxygen concentrators.
Flow meters	
H537-1	Stand-alone low flow meter kit For use with EverFlo oxygen concentrators.
H644	Cabinet integrated pediatric flow meter kit For use with 600, 605, and M600 (non-OPI units only) oxygen concentrators.
528	Locking flow meter kit For use with 600, 605, M600, and M605 oxygen concentrators.
1082784	5 LPM flow meter replacement kit For use with 600, 605, M600, and M605 oxygen concentrators.

Part number	
Gaskets	
1007715	Gasket replacement kit For use with all Millenniums.
1007716	Gasket replacement tool kit For use with 600, 605, M600, and M605 oxygen concentrators.
Humidifier Accessories	
1083699	Humidifier pouch
1101172	SimplyGo humidifier connecting tube kit
1101602	SimplyGo humidifier connecting tube kit and pouch
Miscellaneous	
1118724	User Manual
H649	Caster brake kit For use with EverFlo oxygen concentrators.
920MPC	Charger stand Power cord not included. For use with 920M.
100-0339-00	External spring, single or 30 per case For use with 950 oximeter.
100-0338-00	Gripping insert For use with 950 oximeter.
1017984	Shoulder strap For use with GoLox.
1062966	Lanyard, black, 20" For use with GO2 oximeter
H654	Air inlet retainer kits, 10 per pk. For use with 600 and 605 oxygen concentrators.
1001684	Flow signal conditioner kit
1019440	Reusable wristbands, 10 per pk. For use with WristOx.
1019438	Reusable wristband extenders, 5", 10 per pk. For use with WristOx.
1015581	Silicone grommets, 10 per pk. For use with M10600, and M10605 oxygen concentrators.
Oxygen attachments	
7314	Tubing adapters, 22mm x 22mm, 6 per pk.
312010	O ₂ enrichment attachments, 10 per pk.
312710	O ₂ enrichment attachment
1039642	Humidifier connector tube

Part number	
Oxygen valves	
302418	Pressure valve
302423	Pressure valves, 5 per pk.
1000348	Oxygen valve
Power Supplies	
1116819	SimplyGo Mini DC Power Supply - U.S.
1116818	SimplyGo Mini AC power supply - U.S.
1082664	SimplyGo Mini AC power cord - U.S.
1116830	SimplyGo Mini External Battery Charger - U.S.
1116816	Standard Lithium Ion battery
1116817	Extended Lithium Ion battery
1082662	Lithium Ion battery
1116830	SimplyGo Mini External Battery Charger - U.S.
1083692	SimplyGo DC power supply
1082661	SimplyGo AC power supply
1083693	SimplyGo Airline power cord
1082664	SimplyGo Mini AC power cord - U.S.
1113902	SimplyGo Mini Airline cable
900-102	Rechargeable li ion battery For use with the EverGo oxygen concentrator.
900-103	Battery charger/recalibrator For use with the EverGo oxygen concentrator. For use with the SimplyGo and EverFlo.
900-121	Airline connector For use with the EverGo oxygen concentrator.
1007355	920MP NiMH rechargeable battery
900-104	AC power cord For use with the EverGo oxygen concentrator.
900-101	DC power cord For use with the EverGo (serial number < 109500)
900-120	DC power cord For use with the EverGo (serial number >= 109500)
H617	Power cord replacement kit, 60Hz For use with all Millenniums.
900-105	AC power supply For use with the EverGo oxygen concentrator.



Oxygen accessories (continued)

Part number

Reusable sensors

1037296 Soft sensor, small, 1M

1037297 Soft sensor, medium, 1M

1037298 Soft sensor, large, 1M

Sensor cables

1019444 Sensor adapter cable. For use with WristOx.

Sensor tape

940 Double-sided tape for 934 sensor, 100 per pk.

941 Double sided tape for 935 and 953, 100 per pk.

942 Overwrap tape strips, 100 per pk.

947 Posey wrap for 935 and 953 sensor, 12 per pk. (pictured)

954A Adult FlexiWrap tape strips for 953 sensor, 25 per pk.

954N Pediatric small FlexiWrap tape strips for 935 sensor, 25 per pk.

Sensors

953 Adult sensor, 44 lbs. or more
Includes Flex sensor, 38mm, straight cable,
25 Hydrogel tape strips, and Posey wraps,
1 meter cable.

953-11 Adult sensor, 3 meter cable, 8'

934 Infant sensor, large toe, 4 to 44 lbs.
Includes Flex sensor, 34mm,
with FlexiWrap tape strips.

935 Pediatric small sensor, outstep of foot, 4 lbs. or less
Includes Flex Sensor, 38mm, right angle cable,
with 25 FlexiWrap tape strips and 12 Posey wraps.

936 Adult finger clip sensor, index finger, 66 lbs. or more

936P Pediatric finger clip sensor, index finger,
18 lbs. or more

948 Adult ear clip sensor, lobe of ear, 88 lbs. or more

951 Reflectance sensor middle forehead, 66 lbs. or more

1019522 Finger clip sensor, adult. For use with WristOx.

Ultrafill Cylinders

1113356 MB08 3000 PSI Cylinder with OCD

1113355 MC13 3000 PSI Cylinder with OCD

1113354 MD22 3000 PSI Cylinder with OCD

1113353 ME36 3000 PSI Cylinder with OCD

1116898 MB08 Cylinder with OCD, CANADA

1116899 MC13 Cylinder with OCD, CANADA

1116900 MD22 Cylinder with OCD, CANADA

1116901 ME36 Cylinder with OCD, CANADA

**Accessories are not separately billable to Medicare, except when used with patient-owned equipment. Other payers may reimburse separately for these items. Check with the specific payer.



Home ventilation

Lightweight, versatile, and easy to use, Philips Respironics ventilation technologies are designed for use in the home, hospital and alternative care sites and provide invasive and noninvasive ventilator support for a wide range of adult and pediatric patients.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.



BiPAP AVAPS

Aligned with System One sleep therapy technology, the BiPAP AVAPS noninvasive ventilator provides the comfort of pressure ventilation and the consistent efficacy of assured tidal volume. Adapting automatically to changing patient needs, BiPAP AVAPS incorporates advanced technologies that help simplify patient care.

Part number		HCPCS
DS1160S	BiPAP AVAPS	E0471
DS1160HS	BiPAP AVAPS with humidifier	E0471, E0562
DS1160TS	BiPAP AVAPS, Heated Tube humidifier	E0471, E0562, A4604
CA1161X (Canada)	BiPAP AVAPS	
CA1161XTS (Canada)	BiPAP AVAPS Heated Tube core package	

BiPAP S/T

The BiPAP S/T noninvasive ventilator offers a straightforward user interface and integrated heated humidification for easy setup and management. Digital Auto-Trak automatically adjusts variable trigger cycle thresholds enabling optimum patient/ventilator synchrony while the optional System One Heated Tube humidification allows humidity levels to be tailored to meet patient needs.

Part number		HCPCS
DS1060S	BiPAP S/T	E0471
DS1060HS	BiPAP S/T with humidifier	E0471, E0562
DS1060TS	BiPAP S/T with Heated Tube humidifier	E0471, E0562, A4604

Not available in Canada

Accessories (S/T and AVAPs devices)

Humidification		
DS6H	System One 60 Series humidifier	E0562
IN6H (Canada)		
DS6T	System One 60 Series Heated Tube humidifier (includes 80 watt power supply and Heated Tube)	E0562, A4604
IN6T (Canada)		
1063785	Water chamber (includes both 50 and 60 Series lids) replacement only	A7046
6T-CONVKIT	Heated Tube humidifier upgrade kit	A4604 – Heated Tubing

Tubing

SYS1HT15	System One Heated Tube	A4604
1074400	System One performance tubing (15mm)	A7037
1032907	System One performance tubing (22mm)	A7037

Filters

1035443	Pollen (foam) filter, reusable (1 pk)	A7039
1029330	Pollen (foam) filter, reusable (2 pk)	A7039
1063096	Ultra-fine filter, disposable (1 pk)	A7038
1035442	Ultra-fine filter, disposable (2 pk)	A7038
1029331	Ultra-fine filter, disposable (6 pk)	A7038

Miscellaneous

1099592	Device side panel	
1063857	Carrying case	

To learn more visit www.philips.us/bipapst





Trilogy100 and Trilogy200

Market-leading versatility

Technologies that are designed to automatically adjust to changing patient needs

Lightweight, versatile and easy to use, Trilogy ventilators are designed for use in the home, hospital and alternative care sites and provide invasive and noninvasive ventilator support for a wide range of pediatric and adult patients. They combine Auto-Trak and AVAPS technologies with advanced flow triggering and leak compensation features that allow for a more sensitive delivery of therapy. Trilogy also incorporates mouthpiece ventilation (MPV). A 'kiss' trigger detects when a patient engages and disengages from the mouthpiece to deliver on-demand ventilation. It combines with an optional MPV support system to enhance patient ease of use.

Key benefits

- Mouthpiece ventilation (MPV)
- Pressure support up to 40 cm H₂O
- AVAP-AE
- Adjustable AVAPS
- Auto-Trak and Sensitive Auto-Trak
- On-screen waveforms



Part number		HCPCS
1054260B	Trilogy100 ventilator	E0465 or
CA1054096B	Trilogy100 ventilator, Canada	E0466
1040005B	Trilogy200 ventilator	E0465 or
CA1032800B	Trilogy200 ventilator, Canada	E0466

To learn more visit www.philips.com/trilogy100

Trilogy100 and Trilogy200



MPV circuit support system

An interlocking support system combines with Trilogy's breathing circuit to make access easy for patients. The support system can be attached to various surfaces in many configurations with the use of a clamp designed for use on flat or rounded surfaces. The interlocking support circuit can be adjusted to enhance comfort and accessibility.



Wedge attached

Wedge removed

Part number		HCPSC
Remote alarm		
1118941	Remote Alarm Includes: 50ft RJ9 cable, power supply, 9v battery and manual	
1120633	Power Supply US	
1114933	Remote Alarm Cable, RJ9, 50ft (15 m)	
1114934	Remote Alarm, Cable, RJ9, 200ft (61 m)	
34046	Remote Alarm Stand	
Power cords		
1044760	Power cord, C7-end, RA, FemEnd, 8ft, USA/Can	
1069219	Power cord, C7-end, RA, FedEnd, 14 ft, USA/Can	
Cables		
1045290	Nurse call adaptor cable	
1069979	¼" to 2.5mm NC nurse call cable	
1069980	¼" to 2.5mm NO nurse call cable	
1080247	RJ9 to ¼" nurse call cable alarm closed	
1080248	RJ9 to ¼" nurse call cable alarm open	
1080249	Nurse call RJ9 cable open ends	
1080250	RJ9 to 2.5mm 3" female adapter cable	
1047295	External battery cable	A4612
1067424	Automotive DC power adapter	A4612
Porting blocks		
1054670	Trilogy100, exhalation port block, active	
1040372	Trilogy100, exhalation port block, passive	
1040370	Trilogy200, universal porting block	
Reusable circuits		
1052064	Trilogy active w/PAP exhalation circuit, adult	A4618
1052065	Trilogy active w/PAP exhalation circuit, ped	A4618
1052066	Trilogy passive exhalation circuit, adult	A4618
1052127	Trilogy passive exhalation circuit, pediatric	A4618
1068837	Trilogy200 active w/Flow sensor circuit, adult	A4618
1068838	Trilogy200 active w/Flow sensor circuit, pediatric	A4618
Adult circuits		
1069467	Adult passive circuit with fixed exhalation port, 10 pk	A4618
1069468	Adult passive circuit with fixed exhalation port, 20 pk	A4618
1073228	Adult passive circuit with swivel design exhalation port, 10 pk	A4618
1073236	Adult passive circuit with swivel design exhalation port, 20 pk	A4618
1069443	Adult passive circuit with water trap, 10 pk	A4618
1069444	Adult passive circuit with water trap, 20 pk	A4618
1073227	Adult passive circuit, heated wire, 10 pk	A4618
1073235	Adult passive circuit, heated wire, 20 pk	A4618
1073223	Adult passive circuit, no water trap, 10 pk	A4618
1073232	Adult passive circuit, no water trap, 20 pk	A4618
1073222	Adult active circuit, heated wire, 10 pk	A4618
1073231	Adult active circuit, heated wire, 20 pk	A4618
1102866	MPV circuit, 10 pk	A4618
Pediatric circuits		
1073224	Pediatric passive circuit with water trap, 10 pk	A4618
1073233	Pediatric passive circuit with water trap, 20 pk	A4618
1073225	Pediatric passive circuit, heated wire, 10 pk	A4618
1073234	Pediatric passive circuit, heated wire, 20 pk	A4618
1073220	Pediatric active circuit with water trap, 10 pk	A4618
1073229	Pediatric active circuit with water trap, 20 pk	A4618
1073221	Pediatric active circuit, heated wire, 10 pk	A4618
1073230	Pediatric active circuit, heated wire, 20 pk	A4618

Part number		HCPSC
Circuit adapters		
1073902	Flexible trach adapter with 22 mm cuff (connects to Whisper Swivel II and swivel and fixed disposable exhalation ports), 10 pk	
1073904	Flexible trach adapter with 22 mm cuff (connects to Whisper Swivel II and swivel and fixed disposable exhalation ports), 20 pk	
1073863	Flexible trach adapter with 15 mm cuff (connects to the flow sensor or the active exhalation valve), 10 pk	
1073903	Flexible trach adapter with 15 mm cuff (connects to the flow sensor or the active exhalation valve), 20 pk	
1065659	Active exhalation valve with PAP assembly	
1061398	Trilogy 200, universal active PAP tube adapter	
1050408	Trilogy 200, proximal flow sensor	
1050534	Trilogy 200 Proximal flow sensor, 10 pk	
332113	Whisper Swivel II	
312149	Disposable fixed exhalation port, 10 pk	
1084428	Disposable swivel passive exhalation port, 10 pk	
1084429	Disposable swivel passive exhalation port, 20 pk	
1065775	Filtered exhalation port, 10pk	
Carry cases		
1033709	Trilogy travel carrying case	
1040420	Trilogy carrying case, in-use	
Battery packs		
1072182	Detachable battery charger	A4613
1043570	Trilogy detachable battery pack, dom/intl	A4611
Filters		
1029330	Trilogy filter, 2 pk	
C06417	Filter, bacteria, 22 mm, single use	
C06418	Filter, bacteria, 22 mm, single use, 50 pk	
Stands		
1047410	Trilogy small products roll stand	
1094043	Trilogy circuit support arm	
1074004	RP-mount plate, Trilogy roll stand	
1074137	RP screw, MNT plate, Trilogy roll stand, 4 pk	
Humidifiers		
22135	Humidifier, Fisher & Paykel, MR410	E0550
1074174	Fisher & Paykel airway temperature probe	
1074175	Fisher & Paykel heated hose humidifier adapter	
Miscellaneous		
1102862	MPV support system	
1118840	MPV Replacement Straw Kit (10pk)	
1040390	Trilogy, O2 quick-disconnect insert, 2 pk	
1021671	Test lung	
1045587	Inlet airpath assembly	

Under Medicare, all parts and accessories are included in the monthly rental of the device and are not separately billable. Supplies may be separately billable for patient-owned ventilators. Consult the payer policy manual for specific details.

Trilogy accessories are separately billable only when the ventilator is patient owned.

Airway clearance

For use by caregivers and patients, in the hospital or in the home, CoughAssist T70 offers a gentle and more comfortable alternative to traditional suctioning methods.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.



Clearly a better solution

CoughAssist T70 is a noninvasive therapy that consistently removes secretions in patients with an ineffective ability to cough. It clears secretions by gradually applying a positive pressure (insufflation) to the airway and then rapidly shifting to negative air pressure (exsufflation). This in-exsufflation shift in pressure creates a high expiratory flow that simulates a deep, natural cough. Instead of introducing a suction catheter into the airway, air is delivered noninvasively through a facemask, a mouthpiece, or through an adapter that allows the device to function with an endotracheal or tracheostomy tube.

Key benefits

- Digital display and push button controls
- AC, DC and optional battery power options
- Cough-Trak automatic sensitivity
- DirectView data management



Part number		HCPCS
1098160	CoughAssist T70	E0482
1098159	CA E70 (Canada)	E0482

Circuits		
1090830	CA 70 series circuit, infant, 6 ft	A7020
1090831	CA 70 series circuit, toddler, 6 ft	A7020
1090832	CA 70 series circuit, adult small, 6 ft	A7020
1090833	CA 70 series circuit, adult medium, 6 ft	A7020
1090834	CA 70 series circuit, adult large, 6 ft	A7020
1090835	CA 70 series circuit, trach, 6 ft	A7020
1090836	CA 70 series circuit, mouthpiece, 6 ft	A7020
1098403	CA 70 series circuit, infant, 9 ft	A7020
1098404	CA 70 series circuit, toddler, 9 ft	A7020
1098405	CA 70 series circuit, adult small, 9 ft	A7020
1098407	CA 70 series circuit, adult medium, 9 ft	A7020
1098408	CA 70 series circuit, adult large, 9 ft	A7020
1098409	CA 70 series circuit, trach, 9 ft	A7020
1098410	CA 70 series circuit, mouthpiece, 9 ft	A7020

Masks	
1046861	CoughAssist, facemask, adult large
1044187	CoughAssist, facemask, adult medium
1044188	CoughAssist, facemask, adult small
1044190	CoughAssist, facemask, infant
1044189	CoughAssist, facemask, toddler

Miscellaneous	
1046860	CoughAssist, disposable filter
1099533	CoughAssist, mouthpiece
1098720	CA 70 series, water trap
1098655	CA 70 series, roll stand
1098656	CA 70 series, roll stand bracket
1059017	CA 70 series, foot pedal
1098884	CA 70 series, carrying case
1051801	SD card and case
1099035	CA 70 series, circuit retainer
1032584	Pollen filter, single pack
1043570	Trilogy detachable battery pack, dom/intl
1099036	CA 70 series, battery cover
1072182	Detachable battery charger
1099427	RP-power cord, 8 ft, C17, NEMA 1-15P
1067424	Trilogy/CoughAssist automotive DC power adapter

Accessories are not separately billable to the Medicare program, with the exception of replacement circuits. Consult your local Medicare DME MAC local coverage determination (LCD) for details. Commercial insurance plans may provide separate coverage for accessories. Consult insurance company policy manual for more information.

To learn more about CoughAssist T70 visit www.philips.us/coughassistt70



Respiratory drug delivery

Respiratory Drug Delivery provides innovative drug delivery technologies for the treatment of respiratory diseases. Our easy-to-use products include compressor nebulizer systems, holding chambers, inspiratory muscle trainer, positive expiratory pressure (PEP) device and peak flow meters.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.



InnoSpire Essence

compressor nebulizer system

Compact, lightweight and economical, the InnoSpire Essence compressor nebulizer system combines with proven SideStream technology to provide fast and efficient aerosol drug delivery.

Part number		HCPCS
1099966	InnoSpire Essence with SideStream disposable nebulizer	E0570
1100312	InnoSpire Essence with SideStream disposable and reusable nebulizers	E0570
1100313	InnoSpire Essence with SideStream disposable and reusable nebulizers, Tucker the Turtle mask, and carrying case	E0570

Accessories		HCPCS
Filters		
1114112	Replacement filter, 50 pk	A7013
1102088	Replacement filter, 4 per pk	A7013
Miscellaneous		
1026568	Premium carrying case	
1026641	Standard carrying case	
1101452	InnoSpire carrying case	

To learn more visit www.philips.us/essence



InnoSpire Elegance

compressor nebulizer system

The InnoSpire Elegance compressor combined with the SideStream nebulizer provides a powerful, high-quality compressor nebulizer system that is durable, easy to use, and easy to maintain. It provides the freedom to move, helping your patients build solid compliance behaviors – in a size that easily fits into their lives.

Part number		HCPCS
1100134	InnoSpire Elegance with SideStream disposable nebulizer	E0570
1099969	InnoSpire Elegance with SideStream disposable and reusable nebulizers	E0570
1102828	InnoSpire Elegance with SideStream Plus	E0570

Accessories

Filters

1114112	Replacement filter, 50 pk	A7013
1102088	Replacement filter, 4 per pk	A7013

Miscellaneous

1026568	Premium carrying case	
1026641	Standard carrying case	
1101452	InnoSpire carrying case	



To learn more visit www.philips.us/elegance

InnoSpire Deluxe

compressor nebulizer system

For use with aerosol medications, InnoSpire Deluxe incorporates built-in conveniences that encourage adherence. It includes a convenient storage compartment to hold extra medicine or nebulizer supplies, the attached power cord tucks neatly into the storage area, and a sturdy handle makes it easy to carry. InnoSpire Deluxe is designed with adherence in mind.



Part number		HCPCS
1110063	InnoSpire Deluxe compressor with SideStream disposable nebulizer (6/cs)	E0570
1114112	Replacement filter, 50 pk	A7013
1102088	Filter, 4 pk	A7013





Sami the Seal

compressor nebulizer system

The Sami the Seal pediatric compressor is shaped like a seal to appeal to children, their parents, and pediatricians, and to support aerosol therapy compliance. The compressor partners with highly efficient SideStream nebulizers to provide fast and friendly treatment to pediatric patients in the home.

Part number		HCPCS
1082456	Sami compressor with SideStream disposable nebulizer and Tucker the Turtle mask	E0570
1082731	Sami compressor with SideStream disposable and reusable nebulizers and Tucker the Turtle mask	E0570
Accessories		
1083128	Sami compressor replacement filter (5 pack)	A7013
1044382	Tucker the Turtle pediatric face mask (for use with SideStream reusable H5860 or SideStream disposable H5800)	



InnoSpire Mini

compressor nebulizer system

InnoSpire Mini allows your patients to nebulize their medications whenever they want and wherever they may be. Compact, lightweight and reliable, it provides independence for people who want to get the most out of life.

An ideal traveling companion, InnoSpire Mini's compact size allows for easy transport while a breath-enhanced reusable nebulizer makes it easy to use. The SideStream Plus reusable nebulizer provides fast, effective treatment and can be used to effectively deliver most commonly prescribed aerosol medications.

Choose the time and place that's right for you

An optional rechargeable lithium ion battery operates the InnoSpire Mini compressor for up to 90 minutes between charges (and can be fully recharged in approximately 6 hours). A 12V DC adapter enables the compressor to be powered in autos, trucks, boats and recreational vehicles. The compressor can also be powered by the traditional 110V AC wall outlet. InnoSpire Mini is compliant with selected FAA environmental RTCA/DO-160-F standards for airborne equipment as well as the international basic safety and essential performance standard IEC 60601-1 third edition and its amendments.

Accessories include

- SideStream Plus high-efficiency, reusable nebulizer
- Optional rechargeable lithium ion battery
- 12V DC (car) adapter
- AC adapter
- Carrying case
- Spare filters



Part number

HCPCS

Compressor nebulizer systems

1109289	InnoSpire Mini compressor nebulizer standard	E0570
1109290	InnoSpire Mini compressor nebulizer with battery	E0570

Accessories/replacement parts

1110726	InnoSpire Mini rechargeable lithium ion battery	
1109441	InnoSpire Mini AC power adapter	
1109455	InnoSpire Mini 12V DC car adapter	
1109488	InnoSpire Mini replacement air filter pack	A7013

Aeroneb Go



Aeroneb Go


The Aeroneb® Go nebulizer is a fast, efficient, simple-to-use nebulizer developed for the millions of patients worldwide who require respiratory therapy in and away from home.

The Aeroneb Go can be used by pediatric as well as adult patients and is ideal for conditions such as asthma and chronic obstructive pulmonary disease (COPD). The Aeroneb Go aerosolizes clinician-prescribed solutions for inhalation that are approved for use with a general purpose nebulizer.

Designed with convenience in mind and to help users maintain normal daily schedules, it is ultra-lightweight, portable, quiet and is hassle-free to set up and operate.

Part number	
AG-AG4200-US	Aeroneb Go Complete System: Contains nebulizer unit, battery Pack, wall adapter, extra mouthpiece, mask elbow, IFU, carrying case
AG-AG4100-US	Aeroneb Go Convenience Unit: Contains nebulizer unit, battery pack, extra mouthpiece, mask elbow, IFU, carrying case

Accessories	
AG-AG4300	Aeroneb Go nebulizer handset
AG-AG4500	Aeroneb Go nebulizer handset + connector cable
AG-AG3010	Cap for nebulizer handset (50 pack)
AG-AG3090	Mask Adapter, Intl. 22M / 22F (50 pack)
AG-AG3040	Mouthpiece (50 pack)
AG-AG4060	Control module + connector cable
AG-AG4070	Connector cable (for control module)
AG-AG3410-US	AC/DC mains adapter
AG-AG2110	Carrying case



Sidestream nebulizers

SideStream disposable nebulizer

The SideStream disposable nebulizer offers fast, efficient, and economical drug delivery. It provides consistent drug delivery time after time and can be used with all commonly prescribed respiratory medications. The unique venturi design offers fast drug delivery and short treatment times. It can be used for two weeks before replacing.

Part number		HCPCS
HS800	SideStream disposable nebulizer, 50 pk	A7003
1028219	SideStream disposable cup, 50 pk	A7004
1036160	SideStream disposable with adult mask, 50 pk	A7003, A7015
1041906	SideStream disposable with pediatric mask, 50 pk	A7003, A7015

SideStream reusable nebulizer

The unique venturi design offers fast drug delivery and short treatment times. It can be used for six months before replacing.

Collis G.G. et al. *Lance* 1990;336:341

Part number		HCPCS
HS860	SideStream reusable nebulizer, 10 pk	A7005
HS860-050	SideStream reusable nebulizer, 50 pk	A7005

SideStream Plus reusable nebulizer

The SideStream Plus breath-enhanced nebulizer features an easy-action inspiratory valve that opens on inspiration to boost medication delivery and closes on exhalation, preserving medication and reducing waste. Less pressure is required to open the inspiratory valve, making treatments easier for pediatric and elderly patients and its proven SideStream five-jet venturi design provides a high respirable fraction and shortens treatment times. It can be replaced every 6 months.

Part number		HCPCS
HS870	SideStream Plus breath-enhanced nebulizer, 10 pk	A7005
1110824	SideStream Plus breath-enhanced nebulizer, 50 pack	A7005



Nebulizer Kits

SideStream replenishment nebulizer kit

These custom nebulizer kits are designed to provide a month replacement of product in each package. Includes two complete SideStream disposable nebulizers and two replacement nebulizer cups.

Part number		HCPCS
1028125	SideStream custom nebulizer kits, 25 pk	A7003 and A7004



Nebulizer Masks

SideStream masks

The SideStream mask is recommended for use with the SideStream high-efficiency nebulizers. The reusable mask is contoured to fit the face. The pediatric mask is recommended for ages five and under; the adult size for ages six and up.

Part number		HCPCS
1025529	SideStream pediatric mask, 50 pk	A7015
1025531	SideStream adult mask, 50 pk	A7015
1052521	SideStream Plus adult mask	A7015



SideStream pediatric eye-safe masks

Sami the Seal and Tucker the Turtle child-friendly masks are designed to appeal to children and are constructed to reduce aerosol deposition in and around the eyes. Used with SideStream nebulizers, Sami and Tucker pediatric masks are made of soft, flexible material for comfort and fit.

Part number		HCPCS
HS880	Sami the Seal pediatric mask (for use with SideStream Plus HS870 nebulizer only)	A7015
1044382	Tucker the Turtle pediatric mask	A7015





OptiChamber Diamond

antistatic valved holding chamber

The “one for all” chamber

In the hospital or at home, its uniquely designed stepped mouthpiece and detachable masks make the OptiChamber Diamond suitable for patients of all ages.

The chamber that fits into your patients' lives

OptiChamber Diamond is smaller than most conventional VHCs and is designed for patient comfort. Its stepped mouthpiece facilitates transfer from pediatric mask to mouthpiece and its flat bottom maintains stability when the VHC is not in use.

Key features

- Intuitive, ergonomic design
- Low resistance expiratory valve
- Flexible pMDI adapter
- Detachable LiteTouch facemasks with soft-seal technology

Part number		HCPCS
1077478	OptiChamber Diamond, 10 pk	A4627
1079823	OptiChamber Diamond with LiteTouch small mask, 10 pk	A4627
1079826	OptiChamber Diamond with LiteTouch medium mask, 10 pk	A4627
1079827	OptiChamber Diamond with LiteTouch large mask, 10 pk	A4627
1098213	OptiChamber Diamond, 10 pk, Canada	
1098214	OptiChamber Diamond with Lite Touch small mask, 10 pk, Canada	
1098215	OptiChamber Diamond with LiteTouch medium mask, 10 pk, Canada	
1098216	OptiChamber Diamond with LiteTouch large mask, 10 pk, Canada	




ProChamber valved holding chamber

ProChamber is the professional model valved holding chamber for the acute respiratory care setting. Its proven clinical performance and great value make ProChamber the choice for the budget-conscious clinician.

Part number		HCPCS
HS2003-050	ProChamber valved holding chamber, 50 pk	A4627



Holding chamber masks



LiteTouch valved holding chamber mask


LiteTouch with soft seal technology that attaches onto the mouthpiece of the OptiChamber Diamond valved holding chamber (VHC). Gently touching the mask to the face creates a cushion-like seal that quickly contours to the face, making it easy to fit, especially for small children. Coupled with the Diamond VHC, LiteTouch's soft-seal technology minimizes leaks and helps to deliver optimum aerosol therapy for asthma and other respiratory illnesses. It is not manufactured with BPA as declared by our suppliers.

Part number

1082714 LiteTouch mask small, 10 pk

1082715 LiteTouch mask medium, 10 pk

1082716 LiteTouch mask, large, 10 pk



Silicon masks

Our detachable valved holding chamber masks for the ProChamber provide an effective seal and conform to fit different facial geometries. Built-in reinforcements help prevent the mask from collapsing and assist in maintaining a seal, especially when your pediatric patients are not cooperative. Available in small (0 – 18 months), medium (1 – 6 years), and large (6 years +) sizes.

Part number

HS81111-010 OptiChamber small face mask, 10 pk

HS81211-010 OptiChamber medium face mask, 10 pk

HS81311-010 OptiChamber large face mask, 10 pk

Chambers, and asthma kits

AsthmaPACK

The AsthmaPACK adult and children's personal asthma care kit for patient adherence implements three key elements of the NAEPP asthma treatment guidelines with one prescription-peak flow monitoring, MDI drug delivery, and asthma education. One prescription will give your patients a comprehensive, clinician-directed, patient-managed asthma plan with the tools needed to better understand their asthma, monitor their conditions objectively, and take their prescribed medications correctly and consistently.

Part number

1099014 Adult asthma care kit, 10 pk

1099015 Children's asthma care kit, 10 pk

Kits include a peak flow meter, spacer, instructional video, and educational self-help booklet.



AsthmaPACK children's asthma care kit

AsthmaPACK children's asthma care kit makes learning to manage asthma child's play. This just-for-kids kit makes learning fun with engaging cartoon characters and age-appropriate devices. An exciting feature of this kit is the interactive bilingual (English/Spanish) coloring and activity storybook, "Tucker Tackles Asthma."

Part number

1099014 Adult asthma care kit, 10 pk

1099015 Children's asthma care kit, 10 pk

Kits include a peak flow meter, spacer, instructional video, and educational self-help booklet.



Peak flow meters



AsthmaCheck

AsthmaCheck's easy-grip handle directs hands away from the airflow path and is designed for accurate measurements. Like all of our peak flow meters, it has a built-in three-zone management system and an easy-to-read scale which is suitable to all patients (60 – 810 l/min).

Part number	HCPCS
HS740-012 AsthmaCheck, 12 pk	A4614



Personal Best peak flow meter

Wherever Personal Best goes – in the hospital, office, or home – it sets a new standard for convenience, compliance, and confidence. And Personal Best goes anywhere...easily. A lightweight, self-contained, portable design puts peak flow monitoring at your fingertips – or your patients' – all day long. Plus, a rugged, sterilizable construction helps it stand up to the roughest use day after day. Personal Best encourages compliance with its comfortable oval mouthpiece and easy-to-use built-in handle. And it improves patient communication with its high-legibility scale and NAEPP-endorsed three-zone management system. Available in both full range (60 – 810 l/min) and low range (50 – 390 l/min).

Part number	HCPCS
HS755-012 Personal Best full range peak flow meter, 12 pk	A4614
HS756-012 Personal Best low range peak flow meter, 12 pk	A4614

Mouthpieces

Disposable and one-way valve safety mouthpieces

Mouthpieces help support your goals in reducing cross-contamination during multi-patient screening with Personal Best, and AsthmaCheck peak flow meters. The one-way valve prevents patients from inadvertently breathing in air through their meters during PEFr testing.

Part number	HCPCS
HS715-100 Disposable pediatric plastic mouthpieces for all Respironics' peak flow meters, 100 pk	
HS757-100 Disposable mouthpieces for Personal Best, AsthmaMentor, and AsthmaCheck, 100 pk	
HS714-200 One-way valve safety mouthpieces for Personal Best, and AsthmaCheck, 200 pk	



Inspiratory muscle trainers and airway clearance

Threshold IMT

Provides consistent and specific pressure for inspiratory muscle strength and endurance training, regardless of how quickly or slowly patients breathe. This device incorporates a flow-independent, one-way valve to provide consistent resistance and features an adjustable specific pressure setting (in cm H₂O) to be set by a healthcare professional. When patients inhale through Threshold IMT, a spring-loaded valve provides resistance that exercises respiratory muscles through conditioning.

Part number	HCPCS
HS730-010 Threshold IMT, 10 pk	E1399

Threshold PEP

Threshold PEP is used for airway clearance, bronchial hygiene, or as an alternative to chest physical therapy. Threshold PEP incorporates a flow-independent, one-way valve to ensure consistent resistance with adjustable specific pressure settings (in cm H₂O) which are set by the healthcare professional. When patients exhale through Threshold PEP, the resistive load creates positive pressure that helps open the airways and allows mucus to be expelled during "huff" coughing (forced expiratory technique).

Part number	HCPCS
HS735-010 Threshold PEP, 10 pk	E1399



Sleep diagnostics

Whether you're considering new PSG systems, home sleep testing systems, titration systems, sensors, integrated scoring solutions or other diagnostic tools, we can help you with our full line of diagnostic products, services and solutions. When choosing our diagnostic hardware, we include Sleepware G3 software, unlimited software licensing, a lifetime of core software updates, a starter sensor kit and more to help support and accommodate your business needs.



Somnolyzer

now integrated with Sleepware G3



Somnolyzer from Philips Respironics is a clinically-validated, computer-assisted sleep scoring system designed to help simplify and improve sleep study scoring productivity, accuracy, and consistency. By reducing the time sleep technicians spend on studies – typically by half – they are free to work on other important lab activities.

More details about Somnolyzer include:

- **Somnolyzer increases lab efficiency by reducing the time spent on manual scoring**
Some technicians may find that they only need as little as 10-20 minutes for a structured, high-quality expert review. Somnolyzer can analyze both PSG and HST recordings.
- **Integration with Sleepware enables real-time scoring**
During an in-lab PSG study, Somnolyzer will provide real-time, high quality staging and event scoring. As Somnolyzer applies staging and scoring, the night tech can easily track an ongoing tally of sleep efficiency and sleep time as well as event summaries such as RDI, AHI, and leg movements.
- **Simple to interpret confidence trend or "traffic Light" system**
Sleepware with Somnolyzer uses straightforward green, yellow, and red indicators to show where it has high, medium or low confidence in its' scoring, allowing the expert reviewer to focus their expertise on where it is needed the most.
- **Same scoring rules as human expert to allow for validity**
Somnolyzer recognizes sleep characteristics — spindles, eye movements and delta waves — and applies the same decision process to staging and event scoring as a human expert would. The Somnolyzer's scores are indistinguishable from those of a human scorer.
- **Diminishes interscorer variability for increased reliability**
Somnolyzer scoring, along with a structured expert review, diminishes inter-scorer variability.
- **Follows R&K and AASM standards to deliver excellent results**
Patient results are compliant with the latest professional guidelines.

Somnolyzer is fully integrated within Sleepware G3, and is an available option with a per-use license activation.

Description	Part no.
Somnolyzer-G3 Act., Full, 20 uses	1121056
Somnolyzer-G3 Act., Full, 50 uses	1121057
Somnolyzer-G3 Act., Full, 150 uses	1121148
Somnolyzer-G3 Act., Full, 250 uses	1121149
Somnolyzer-G3 Act., Full, 500 uses	1121150
Somnolyzer-G3 Act., Full, 1000 uses	1121151
Somnolyzer-G3 Act., Full, 2000 uses	1121152
Somnolyzer-G3 Act., Resp Only, 20 uses	1121153
Somnolyzer-G3 Act., Resp Only, 50 uses	1121154
Somnolyzer-G3 Act., Resp Only, 150 uses	1121155
Somnolyzer-G3 Act., Resp Only, 250 uses	1121156
Somnolyzer-G3 Act., Resp Only, 500 uses	1121157
Somnolyzer-G3 Act., Resp Only, 1000 uses	1121158
Somnolyzer-G3 Act., Resp Only, 2000 uses	1121159

Alice 6 LDx

Sleep system

The LDx base station can be paired with either the LDxS or LDxN Headbox, so your sleep lab can have the right capability where and when you need it. The LDxS is the perfect front-line headbox for most sleep center needs, with 19 EEG inputs. The LDxN Headbox offers 32 EEG inputs. Both head boxes are easy-to-understand with single keyhole inputs allowing staff to focus more on patients and less on wires and equipment.

Alice 6 LDx base system,

North America

1063314

Includes:

- Alice 6 LDx base station
- AC power cord
- Auto-switching power supply
- Alice 6 communication cable
- Alice 6 user guide
- Headbox shoulder strap (headbox not included)
- Unlimited site-licensing
- 2-year warranty
- Software updates for life

Alice 6 LDx base system,

International

1063315

Includes:

- Alice 6 LDx base station
- AC power cord
- Auto-switching power supply
- Alice 6 communication cable
- Alice 6 user guide
- Headbox shoulder strap (headbox not included)
- Unlimited site-licensing
- 2-year warranty
- Software updates for life

LDxS Headbox

1063552 North America

1063353 International



LDxN Headbox

1063476 North America

1063547 International



To learn more visit www.philips.com/sleepdx



Alice 6 LDE

Sleep system

The Alice 6 LDE diagnostic system provides your busy lab with the base set of channels needed to meet AASM standards while allowing staff to focus more on patients and less on equipment.

LDE Headbox
1063549 North America

LDE Headbox
1063550 International



Alice 6 LDE base system, **North America** **1063311**

Includes:

- Alice 6 LDE base station
- AC power cord
- Auto-switching power supply
- Alice 6 communication cable
- Alice 6 user guide
- Headbox shoulder strap (headbox not included)
- Unlimited site-licensing
- 1-year warranty
- Software updates for life

Alice 6 LDE base system, **International** **1063312**

Includes:

- Alice 6 LDE base station
- AC power cord
- Auto-switching power supply
- Alice 6 communication cable
- Alice 6 user guide
- Headbox shoulder strap (headbox not included)
- Unlimited site-licensing
- 1-year warranty
- Software updates for life

To learn more visit www.philips.com/sleepdx

Computers/monitors



Description	Part no.
LCD computer monitor, 20"	1013805
LCD computer monitor, 24"	1043655
Computer, desktop for Alice systems	1043167
Computer, laptop for Alice systems	1024624
Speakers for Alice PSG	1118288
CPU/laptop microphone for Alice PSG	1025779
Microphone kit, super-high gain	1097360

Sleepware G3 with Somnolyzer software

One software platform that powers both our in-lab PSG systems and portable devices.

Description	Part no.
Sleepware G3, English only	1082462
Sleepware G3, Multi-language	1082463
HL7 inbound only interface	1069193
HL7 outbound only interface	1069194
HL7 bi-directional interface	1067148

Video/audio



Description	Part no.
Video camera (fixed focus, analog)	1016348
Axis Video Server (single port),	1071285
Analog video cable, 50 ft.	1021731
Analog video cable, 100 ft.	1021730
VideoEase CCTV balun	1026576
Modular CCTV power – through balun	1032204
Axis pan-tilt-zoom network camera	1117858
IR projector	1016349
Microphone kit, super high-gain	1097360
Speakers for Alice PSG	1094032

Alice lab system accessories

Description	Part no.
CPAP therapy circuit interface kit, Alice 6	1082653
Headbox communication cable, Alice 6	1074689
Wall mounting kit, Alice 6	1082652
Shoulder strap for Alice 6 headbox	1078921
Alice 5 headbox wall mounting kit	1020266
Alice 5 headbox carrying pouch	1020584
Shoulder strap for Alice 5 headbox	1016343
Headbox communication cable, Alice 5	1024583
Sleepware G3 Getting Started Guide, English	1081398
Sleepware user guide, English	1016318
Sleepware user guide, English/International	1016319

Cables and network



H3150-50, H3150-100 H3151-50, H3151-100 1024583 8158-10 1016352 1016351

1007492 1021726, 1021727 8179 8155

Description	Notes	Part no.
Headbox communication cable, Alice 6		1074689
Connecting cable CAT-5 patch, 50'	Standard network patch cable	H3150-50
Connecting cable CAT-5 patch, 100'	Standard network patch cable	H3150-100
Connecting cable CAT-5 cross-over, 50'	Cable used for direct connections only, (i.e., PC to base station)	H3151-50
Connecting cable CAT-5 cross-over, 100'	Cable used for direct connections only, (i.e., PC to base station)	H3151-100
Patient cable, Alice 5	Connects base station to headbox	1024583
Sub-mini adapter for auxiliary input cable		8158-10
Connecting cable auxiliary input, 10'		1016352
A/B aux input signal splitter		1016351
Reference electrode jumper cable		8179



1012148

Description	Part no.
Network Switch, Gigabit, 16-port	1106684

Miscellaneous accessories




1024582 1016309 1016310 1016344

Description	Part no.
Photo detector	1050795
Hard shell carrying case includes: sensor case	1024582
Sensor case, sensors not included	1016309
Mouse pad	1016310
Power supply, 60W, 6.25 VDC 3-6 AMPS, Alice 6	1019280
Uninterruptible power supply	1091779
Power Cord, US domestic (used with Alice 5 and Alice 6)	1120187
Power Supply, (used with Alice 5 and Alice 6)	1019280
Power Supply and Power cord combined (used with Alice 5 and Alice 6)	10224477L

Manuals

Description	Part no.
Alice 6 hardware setup guide, English	1078160
Alice 6 hardware setup guide, English/International	1078161
Alice 5 hardware setup guide, English	1016329
Alice 5 hardware setup guide, English/International	1016330

Sensors – air flow



Description	Notes	Part no.
Thermistor airflow, adult, 1 channel, 249, Alice 6		1078754
Thermistor airflow, pediatric, 1 channel, 249, Alice 6		1078755
Thermistor airflow, pediatric small, 1 channel, 249, Alice 6		1078756
Thermistor, oral, adult, 249, Alice 6		1078757
Thermistor airflow, adult, 1 channel, 1.5mm, Alice 5		P1274
Thermistor airflow, pediatric, 1 channel, 1.5mm, Alice 5		P1336
Thermistor airflow, pediatric small, 1 channel, 1.5mm, Alice 5		P1337
Thermistor, oral, adult, 1.5mm, Alice 5		P1344
Pressure sensor, North America, Alice 6		1076590
Pressure sensor, International, Alice 6		1076592
PTAF Lite, Alice 5		P1304R
Nasal cannula, adult (Pro-Flow)	60 per pk	P1259
Nasal cannula, adult (Pro-Flow)	10 per pk	P1301
Nasal/oral cannula, adult (Pro-Flow Plus)	30 per pk	P1295
Nasal/oral cannula, adult (Pro-Flow Plus)	5 per pk	P1300
Nasal cannula, pediatric, small (Pro-Flow)	60 per pk	P1267
Nasal cannula, pediatric, small (Pro-Flow)	10 per pk	P1303
Nasal cannula, pediatric (Pro-Flow)	60 per pk	P1263
Nasal cannula, pediatric (Pro-Flow)	10 per pk	P1302

Sensors – effort



Description	Notes	Part no.
zRIP, DuraBelt kit, adult, Alice 6	Kit contains: 2 ea. P/N 1073622 zRIP DuraBelt adult effort belt and 2 ea. P/N P1133 wiresets	1078925
zRIP, DuraBelt kit, pediatric, Alice 6	Kit contains: 2 ea. P/N 1075970 zRIP DuraBelt pediatric effort belt and 2 ea. P/N P1133 wiresets	1078926
zRIP Inductance Effort, Alice 5/LE Adult	Kit contains: 2 ea. P/N 1073622 zRIP DuraBelt adult effort belts, P1736 zRIP Alice 5 driver module, 2 ea. P/N P1133 wiresets	1076333
zRIP Inductance Effort, Alice 5/LE Pediatric	Kit contains: 2 ea. P/N 1075970 zRIP DuraBelt pediatric effort belt, P1736 zRIP Alice 5 driver module, 2 ea. P/N P1133 wiresets	1076334
zRIP effort belt, adult		1073622
zRIP effort belt, pediatric		1075970
zRIP effort belt, large		1075973
zRIP module, Alice 5		P1736R
zRIP wireset		P1133
ezRIP adult kit, Alice 5/LE, 84"	Kit contains: 1 abd ezRIP module, and 1 thx ezRIP module, 2 zRIP DuraBelts, adult	1099515
ezRIP module, abdomen, 1.5, 84", Alice 5/LE		1094110
ezRIP module, thorax, 1.5, 84", Alice 5/LE		1094111
CT2, adult, Alice 5	Piezo	P1586
CT2, pediatric, Alice 5	Piezo	P1587

Sensors – oximetry

Description	Part no.
Nonin Xpod SpO ₂ module for Alice 6	1122672
SpO ₂ Sensor (Finger Clip)	936
SpO ₂ Sensor (Flex)	953
Adult Flexwrap (Disposable for use w/953 sensor, 25/pkg)	954A
Nonin Finger Soft Sensor, Large	1062049
Nonin Finger Soft Sensor, Medium	1062050
Nonin Finger Soft Sensor, Small	1062051
Disposable sensor, Adult, cloth, 24 ct	1070517
Disposable sensor, Ped, cloth, 24 ct	1070518
Disposable sensor, Infant, cloth, 24 ct	1070520

Sensors – body movement



P1694

1009963

1025751

Description	Notes	Part no.
Body position sensor, 249, Alice 6		1078758
Body position sensor, 1.5mm, Alice 5		P1694
PLM sensor, 249, Alice 6		1078752
Actimeter, Alice 5	Connects to the wrist, detects body movement	1009963
Actimeter tape	10 yds per roll	1025751

Sensors – electrodes



1016360

1016361

1075397

1075402

8179

8051

8052

Description	Part no.
ECG, 72", DOM pack of 5 lead wires, domestic	1075232
ECG, 72", DOM pack of 5 lead wires, international	1075233
Adult repositionable ECG electrodes	1016360
Infant ECG electrodes (1 pkg. of 3 electrodes)	1016361
Reference electrode jumper cable, 10 ft.	8179
EEG electrode cream	8051
EEG electrode prep	8052
EEG/EMG/EOG, 48", gold cup	1075294
EEG/EMG/EOG, 72", gold cup	1075397
EEG/EMG/EOG, 96", gold cup	1075398
EEG/EMG/EOG, 120", gold cup	1075399
EMG, 72", snap, 10 per pack	1075400
EMG, 96", snap, 10 per pack	1075401
EMG, 120", snap, 10 per pack	1075402

Sensors – snore



P1716, P1696, 1078753

Description	Notes	Part no.
Snore microphone, 249, Alice 6		1078753
Snore microphone, 1.5mm, Alice 5		P1716
Piezo snore sensor		P1696

OmniLab Advanced +

Multi-mode titration system and accessories

With nine therapy modes, a quiet platform, customizable reports, user-friendly features and value-added benefits, you'll find the OmniLab Advanced + to be an ideal system for satisfying your lab titration needs. Along with OmniLab Direct and OmniLab Connect, the OmniLab Advanced + system provides a cohesive, harmonized solution to help you and your sleep technicians operate titration studies efficiently, easily and intuitively.



Description	Part no.
OmniLab Advanced + with Heated Tube humidifier	1111124
A Series humidifier	1111551
Water chamber (includes both 50 and 60 Series lids)	1063785
System One Heated Tube, Performance Tubing (22mm)	1032907
OmniLab Connect (for use with non Philips Respironics PSG systems)	1111076
OmniLab Direct (lab titration software)	1112033
Link Module	1061644A
Serial Cable DB9 Female to Male 6ft	1037268

LoFlo capnography sensor

LoFlo connects to an Alice PSG system with a single USB connection and reduces the need to manage bulky equipment and calibration settings while patients wait to begin their sleep study.



Description	Part no.
Alice, LoFlo C5 Sidestream module	1097955
Alice LoFlo adapter cable	1098078
Sampling line, 10 pk	3474-00
Alice, LoFlo starter kit, (LoFlo, sampling line 10 pk, Alice LoFlo USB cable, and cannulas)	1101087
Cannula, Pro-Flow, nasal, adult 10 pk	1097345
Cannula, Pro-Flow, nasal, pediatric 10 pk	1097348
Cannula, Pro-Flow, nasal, pediatric, small 10 pk	1097349
Cannula, Pro-Flow, nasal, pediatric 60 pk	1097939
Cannula, Pro-Flow, nasal, pediatric, small 60 pk	1097940




Alice NightOne

Home sleep testing

Keep it simple for patients

- “Auto-on” feature: recording starts as soon as device is buckled
- “Smart guide” sequence: easily guides setup and indicates signal quality
- Helpful visual aids: tutorial video; step-by-step diagram



Description	Part no.
Alice NightOne North America Wireless	1113277
Alice NightOne International Wireless	1114801
Alice NightOne International (without Bluetooth)	1114803

Replacement accessories	
Thoracic effort belt, adult	1113310
Thoracic effort belt, large	1113311
Pro-Flow nasal cannula, adult, qty. 60	P1328-60
Carrying case	1117455
Case insert replacement kit	1114520
Quick Start guide, 10 pack	1113329

The Alice NightOne home sleep testing (HST) solution combines ease, comfort and reliability to support your business, clinical and patient needs when it comes to HST. With features and functions designed with patients in mind, Alice NightOne is up to two times more likely to result in an acceptable study than competitive home sleep testing devices.*

Note: Alice NightOne Provider manuals are available in various languages. To order a manual, please contact customer service.

- Key benefits**
- “Auto-on” feature
 - RIP with integrated buckle design
 - Unique “smart guide” patient setup sequence
 - Enhanced Good Study Indicator
 - Body Position sensor built-in
 - Powered by Sleepware G3 software
 - Optional Cloud services
 - Bluetooth connectivity to System One or Dream Station therapy devices to record therapy channels

- Alice NightOne includes:**
- Alice NightOne device
 - Hard shell carrying case
 - Nasal cannula sample pack
 - zRIP effort belt
 - Oximetry module
 - SpO₂ sensor (sensor style selected upon order)
 - USB cable
 - Quick Start guide, user and provider manuals
 - Sleepware G3 software
 - Unlimited software updates and site licensing
 - Two-year device warranty

- Included is a choice of one of the following:**
- 936 SpO₂ Sensor (Finger Clip)
 - 953 SpO₂ Sensor (Flex)
 - 954A Adult Flexwrap (Disposable for use w/953 sensor, 25 pack)
 - 1062049 Reusable Finger Sensor (Soft), Large
 - 1062050 Reusable Finger Sensor (Soft), Medium
 - 1062051 Reusable Finger Sensor (Soft), Small



*Results of customer preference trial; data on file.

**Available in selected countries.

Alice PDx

Home sleep testing

Alice PDx is a portable at-home diagnostic recording device for OSA screening, follow-up, and diagnostic assessment of cardio-respiratory sleep disorders. The system satisfies the AASM/CMS HST device guidelines and enables basic screening to comprehensive diagnostic evaluation. The Good Study Indicator conveys the amount of good quality data collected overnight to help you avoid or reduce the costs associated with retesting.



Description	Part no.
Alice PDx device kit, (Domestic/US)	1043941
Alice PDx device kit, (International)	1043844
Alice PDx device kit, (Japan)	1047990
Alice PDx device kit, (China)	1047992

Alice PDx device kits (batteries, SpO2 probe, and SD card reader, Sleepware also included)

Optional ExG expansion yoke



1040809

Description	Notes	Part no.
ExG yoke, Alice PDx (US and Japan)	Includes 1 ea. of P/N 1020971 and 1016339	1040808
ExG yoke, Alice PDx (International)	Includes 1 ea. of P/N 1020971 and 1016339	1048415
EEG/EMG/EOG leads 72 in., qty 10		1016339
EEG/EMG/EOG leads 48 in., qty 10		1020971

Expansion cable

Description	Notes	Part no.
Alice PDx PLM Sensor	Records movement, uses ExG/Actimeter port)	1069622

Optional ECG expansion yoke



1040809, 1040810,
1048916



1041199



1041200



1016360

Description	Notes	Part no.
ECG yoke, Alice PDx (US)	Includes 1 ea. of P/N 1040809 and 5 ea. of P/N 1016360 and Sleepware analysis	1040809
ECG yoke, Alice PDx (International)	Includes 1 ea. of P/N 1041200 and 5 ea. of P/N 1016360 and Sleepware analysis	1040810
ECG yoke, Alice PDx (Japan)	Includes 1 ea. of P/N 1041200 and 5 ea. of P/N 1016360 and Sleepware analysis	1048916
ECG lead set 5-wire, 30 in. (76 cm), Domestic/US		1041199
ECG lead set 5-wire, 30 in. (76 cm), International		1041200
Adult ECG electrodes (pkg. of 5 electrodes)		1016360

Alice PDX Sensors, Accessories and replacement parts



Description	Notes	Part no.
zRIP DuraBelt adult sensor kit, Alice PDX		1075969
Oral thermistor, Alice PDX		P1379
SpO ₂ extension cable		116888
SpO ₂ sensor (finger clip)		936
SpO ₂ sensor (flex)		953
Adult flexwrap (disposable for use w/953 sensor, 25/pkg)		954A
zRIP DuraBelt effort sensor, adult single belt		1073622
zRIP wireset, Alice PDX		P1170
Alice PDX Thermistor, Airflow, Nasal Oral, Adult, Sensor		P1388
Pro-Flow nasal cannula, adult, 16", 60 per case		P1328-60
Pro-Flow Plus nasal/oral cannula, adult, 16", 30 per case		P1343
CPAP titration kit, Alice PDX		P1391
Alice PDX LM sensor	Includes 1 ea. Alice PDX LM sensor, 1 Velcro strap PLM, small (20cm/18"), and 1 Velcro strap PLM, large (41cm/16")	1069622
Alice PDX holster		1053280
Alice PDX carry case		1043937
Secure digital card 1GB (2 SD cards per package)		1069860
Lanyard (replacement support lanyard)		1053194
SD card reader		1047300
Nonin finger soft sensor, large		1062049
Nonin finger soft sensor, medium		1062050
Nonin finger soft sensor, small		1062051
Surgical tape (box of 12 rolls of 10 yards ea.)		1065281

Sleepware G3 with Somnolyzer software

One software platform that powers both our in-lab PSG systems and portable devices.

Description	Part no.
Sleepware G3, English only	1082462
Sleepware G3, Multi-language	1082463
HL7 inbound only interface	1069193
HL7 outbound only interface	1069194
HL7 bi-directional interface	1067148

Communication and therapy cables

Description	Part no.
Alice PDX sleep therapy cable	1040805
Alice PDX SleepLink cable	1040806
Alice PDX MobileLink cable	1040807
SleepLink card	1059963

Manuals – professional

Description	Part no.
Alice PDX professional guide, (Domestic/US)	1059953
Alice PDX professional guide, (English International)	1050119
Alice PDX professional guide, (French)	1055341
Alice PDX professional guide, (Italian)	1055342
Alice PDX professional guide, (German)	1055343
Alice PDX professional guide, (Spanish)	1055344
Alice PDX professional guide, (Brazilian Portuguese)	1055345
Alice PDX professional guide, (Japanese)	1055346

Alice PDX Professional and User/Patient manuals are available in various languages, please call customer service to order.

Manuals – user/patient

Description	Part no.
Alice PDX user/patient guide, (Domestic/US)	1059952
Alice PDX user/patient guide, (English International)	1050121
Alice PDX user/patient guide, (French)	1055347
Alice PDX user/patient guide, (Italian)	1055348
Alice PDX user/patient guide, (German)	1055349
Alice PDX user/patient guide, (Spanish)	1055350
Alice PDX user/patient guide, (Brazilian Portuguese)	1055351
Alice PDX user/patient guide, (Japanese)	1055352

Stardust 2 Sensors



1094097

Description	Notes	Part no.
Effort ezRIP module, thorax, harwin, 26"		1094097



1007199



612054



1101838

Description	Notes	Part no.
Oximetry Oximeter finger clip sensor	Replaces Y sensor	1007199
Flexible finger oximetry Y probe kit	Includes adult sensor tape	612054
Butterfly wrap, adult	10 per pk	1101838



P1396



P1340



P1328-60



P1343

Description	Notes	Part no.
Airflow Pro-Flow adult nasal cannula, 16"	60 per pk	P1328-60
Pro-Flow Plus adult nasal/oral cannula, 16"	30 per pk	P1343
Stardust PTAF Sensor Kit, description	Includes 5 adult nasal cannulas, 16"	P1396
Stardust PTAF Sensor		P1340

Miscellaneous



1008312



1003510

Description	Part no.
Pouch	1008312
Respironics screwdriver	1003510


Manuals

Description	Part no.
Host software manual v2.0, Domestic	1026567
Host software manual v2.0, International English	1026566
Host software manual v2.0, French	1018899
Host software manual v2.0, German	1026564
Host software manual v2.0, Spanish	1026565



Pro-Tech sensors

As a leader in the development, manufacturing, and supply of sleep diagnostic sensors, we offer a full line of high-performance, quality Pro-Tech sensors to help you maintain low cost of ownership while obtaining the data you need. Our sensor product line includes durable and reliable sensors that detect respiratory airflow, respiratory effort, snoring, periodic limb movements, and sleep position.



Pro-Flow and Pro-Flow Plus

Nasal cannulas

- Designed to enhance nasal and oral signal transfer
- Soft and narrow nasal prongs help reduce nasal resistance and increase patient comfort
- Incorporate narrow kink-resistant tubing to help facilitate excellent response to high-frequency signals such as snoring
- Hydrophobic, anti-microbial inline filter
- Pediatric and adult sizes

Description	Part number
Pro-Flow, nasal, adult, 7', (60 per case), male luer lok	P1259
Pro-Flow, nasal, adult, 7', (10 pack), male luer lok	P1301
Pro-Flow, nasal, adult, 8.25", (60 per case), male luer lok	P1306
Pro-Flow, nasal, adult, 8.25" (30 per case)	P1308
Pro-Flow, nasal, adult, 16", (60 per case), male luer lok	P1328-60
Pro-Flow, nasal, adult, 16", (5 pack), male luer lok	P1390
Pro-Flow, nasal, adult, 7', (10 pack), no filter, male luer lok	1091077
Pro-Flow, nasal, adult, 7', (60 per case), no filter, male luer lok	1091078
Pro-Flow, nasal, pediatric, 7', (60 per case), male luer lok	P1263
Pro-Flow, nasal, pediatric, 7', (10 pack), male luer lok	P1302
Pro-Flow, nasal, pediatric, 16", (60 per case), male luer lok	P1381
Pro-Flow, nasal, Cannula, 8.25" (30 per case), no filter	P1334
Pro-Flow, nasal, pediatric small, 7', (60 per case), male luer lok	P1267
Pro-Flow, nasal, pediatric small, 7', (10 pack), male luer lok	P1303

Pro-Flow Plus


Nasal/oral cannulas

- Addresses monitoring of mouth-breathing patients
- Oral prong detects pressure changes in the mouth
- Pediatric and adult sizes

Description	Part number
Pro-Flow Plus, nasal/oral, adult, 7', (30 per case), male luer lok	P1295
Pro-Flow Plus, nasal/oral, adult, 7', (5 pack), male luer lok	P1300
Pro-Flow Plus, nasal/oral, adult, 8.25", (30 per case), male luer lok	P1308
Pro-Flow Plus, nasal/oral, adult, 8.25", (30 per case), no filter, male luer lok	P1334
Pro-Flow Plus, nasal/oral, adult, 1.1M, (10 per case), male luer lok	P1319-10
Pro-Flow Plus, nasal/oral, adult, 16", (30 per case), male luer lok	P1343
Pro-Flow Plus, nasal/oral, pediatric, 7', (30 per case), male luer lok	P1354



Description	Part number
Pro-Flow MultiPurpose nasal cannulas	
Pro-Flow MultiPurpose, nasal, adult, 7', (30 per case), male luer lok	P1350
Pro-Flow MultiPurpose, nasal, adult, 7', (5 pack), male luer lok	P1368
Pro-Flow MultiPurpose, nasal, pediatric, 7', (30 per case), male luer lok	P1346
Pro-Flow MultiPurpose, nasal, pediatric small, 7', (30 per case), male luer lok	P1348
Pro-Flow Plus MultiPurpose cannulas	
Pro-Flow Plus MultiPurpose, nasal/oral, adult, 7', (30 per case), male luer lok	P1373
Pro-Flow Plus MultiPurpose, nasal/oral, pediatric, 7', (30 per case), male luer lok	P1366
Cannula accessories	
Cannula adapter kit, female luer thread style on 1/8 inch tubing	P1324
Double-port CPAP titration kit Kit contains: 12 cannula titration filters, 2 male luer loks, 2 circuit connection ports, small bore connecting tubing	P1327
CPAP titration kit, Alice PDx, connects CPAP circuit to PDx via male luer lok	P1391
Cannula or titration filter, (50 pack), male and female luer lok	P1107
Cannula or titration filter, (10 pack), male and female luer lok	P1109
Pro-Wrap, 100 Pro-Wrap per roll, no active sensor	P1162
Pro-Flow sample pack, includes P1259 (qty: 2) and P1295 (qty: 2)	P1258



Pro-Flow MultiPurpose and Pro-Flow Plus MultiPurpose cannulas

- Designed to concurrently monitor airflow and end-tidal CO₂
- Provides ability to monitor three parameters (airflow, snore, and EtCO₂) with one patient interface
- Can be used with most pressure transducers and continuous sampling capnographs
- Hydrophobic, anti-microbial inline filter
- Pediatric and adult sizes

Thermal airflow



Thermocouples and thermistors

- Compatible with most major PSG systems
- Ability to monitor oral/nasal separately and/or together with adult and pediatric sizes and disposable thermocouples
- Insulated tips
- Six-month warranty

Description	Part number
Reusable thermocouple (cannula style)	
Thermocouple 6" single element, adult, 1.5mm connectors	P1221
Thermocouple nasal/oral airflow sensor, 1 channel, adult, 1.5mm connectors	P1222
Thermocouple nasal/oral airflow sensor, 2 channel, adult, 1.5mm connectors	P1224
Thermocouple airflow sensor, 1 channel, adult, binder connector	P1757
Thermocouple airflow sensor, 1 channel, adult, binder connector	P1758
Thermocouple nasal/oral airflow sensor, 1 channel, pediatric, 1.5mm connectors	P1223
Thermocouple airflow sensor, 1 channel, pediatric, binder connector	P1759
Thermocouple airflow sensor, 1 channel, pediatric, binder connector	P1760
Thermocouple nasal airflow sensor, 1 channel, pediatric small, nasal bridgemount, 1.5mm connectors	P1214
Thermocouple airflow sensor, 1 channel, pediatric small, binder connector	P1761
Thermocouple airflow sensor, 1 channel, pediatric small, binder connector	P1762
Disposable thermocouple	
Reusable wireset	P1238
Disposable thermocouple nasal/oral airflow sensor, adult, (10 pack), 1.5mm connectors	P1242

Description	Part number
Reusable thermistor	
Alice 3	
Thermistor, 1 channel, nasal/oral, adult, 2.5 connector	P1246
Thermistor 1 channel, nasal/oral, pediatric, 2.5mm connector	P1298
Alice 4	
Thermistor, 1 channel, nasal/oral, adult, LEMO triaxial connector	P1286
Alice 5	
Thermistor, 1 channel, nasal/oral, adult, 1.5mm connectors	P1274
Thermistor 1 channel, nasal/oral, pediatric, 1.5mm connectors	P1336
Thermistor 1 channel, pediatric small, nasal bridgemount, 1.5mm connectors	P1337
Alice 6	
Thermistor, 1 channel, nasal/oral, adult, 249 connector (2-pin keyhole)	1078754
Thermistor, 1 channel, nasal/oral, pediatric, 249 connector (2-pin keyhole)	1078755
Thermistor 1 channel, pediatric small, nasal bridgemount, 249 connector (2-pin keyhole)	1078756
Thermistor, oral, adult, 249 connector (2-pin keyhole)	1078757
Alice PDx	
Thermistor airflow, nasal/oral, adult, Alice PDx, 2.5mm connector	P1388
Thermistor, oral, adult, Alice PDx, 2.5mm connector	P1379
Other	
Thermistor, 1 channel, nasal/oral, adult, Sandman SD20, 220 connector (black)	P1364
Thermistor, 1 channel, nasal/oral, adult, Embletta, 249 connector (2-pin keyhole) (gray)	P1274-249
Thermistor, 1 channel, nasal/oral, pediatric, Embletta, 249 connector (2-pin keyhole) (gray)	P1298-249



Pressure transducers



Thermocouples and thermistors

- Designed to generate accurate respiratory signals
- Includes high-level and low-level outputs, unfiltered or filtered airflow output, snoring output, and remote-control options
- Used in AC or DC configurations and during CPAP titrations
- One-year warranty

Description	Part number
PTAF2	
PTAF2 (kit contains: PTAF2 module, a 9V battery, and one P1258 cannula sample pack)	P1287
PTAF2 cables	
High level cable, 7', 1/8 male connector	P1078
High level cable, 7', 1/8 female connector	P1079
High level cable, 7', 3/32 male connector	P1080
High level cable, 7', 1/4 male connector	P1081
High level cable, 7', 3.5mm connector	P1138
Low level cable, 7', 3/32 male connector	P1077
Low level cable, 7', 1.5mm connector	P1075

Description	Part number
Pressure sensor	
Pressure sensor, domestic, Alice 6	1076590
Pressure sensor, international, Alice 6	1076592



PTAF Lite	
PTAF Lite startup kit	P1289
Kit contains: P1258: cannula sample pack; P1289: PTAF Lite	
PTAF Lite for Alice 5	P1304R

PTAF2 and PTAF Lite accessories	
Cannula titration filter (50 pack)	P1107
Cannula titration filter (10 pack)	P1109
Cannula adapter kit	P1324
CPAP titration kit	P1327
CPAP titration kit Alice PDx	P1391
Cannula sample pack	P1258



The PTAF Lite is a small, lightweight piezo-style pressure transducer used in AC device configurations.

Respiratory inductance plethysmography effort sensors



zRIP respiratory effort sensors and zRIP DuraBelts

- Provide linear breathing movement tracings with no false paradoxing
- Incorporate enhanced woven wire and buckle design
- Offer comprehensive patient coverage
- Facilitate easy initial placement, and easy middle-of-the-night adjustments
- Fully washable
- Available in non-SUM and SUM RT (breath-to-breath summation) modules
- Replaceable batteries
- One-year warranty

Description	Part number
zRIP DuraBelt startup kits	
Alice 5	
zRIP DuraBelt kit, adult, Alice 5 Kit contains: 1 each P1736R zRIP Alice 5, 2 each P1133 zRIP wiresets, 2 each 1073622 zRIP DuraBelts, adult	1076333
zRIP DuraBelt kit, pediatric, Alice 5 Kit contains: 1 each P1736R zRIP driver module, 2 each P1133 zRIP wiresets, 2 each 1075970 zRIP DuraBelts, pediatric	1076334
Alice 6	
zRIP DuraBelt kit, adult, Alice 6 Kit contains: 2 each 1073622 zRIP DuraBelts, adult, 2 each P1133 zRIP wiresets	1078925
zRIP DuraBelt kit, pediatric, Alice 6 Kit contains: 2 each 1075970 zRIP DuraBelts, pediatric, 2 each P1133 zRIP wiresets	1078926
Alice PDx	
zRIP DuraBelt kit, adult, Alice PDx Kit contains: 2 each P1170 wiresets, 2 each 1073622 zRIP DuraBelt, adult	1075969
Other	
zRIP DuraBelt kit, adult Kit contains: 1 each P1732 zRIP driver module, 2 each P1133 zRIP wiresets, 2 each 1073622 zRIP DuraBelt, adult	1076329
zRIP DuraBelt kit, pediatric Kit contains: 1 each P1732 zRIP driver module, 2 each P1133 zRIP wiresets, 2 each 1075970 zRIP DuraBelts, pediatric	1076330
zRIP DuraBelt Sum RT kit, adult Kit contains: 1 each P1753 zRIP Sum RT driver module, 2 each P1133 zRIP wiresets, 2 each 1073622 zRIP DuraBelts, adult	1076336
zRIP DuraBelt Sum RT kit, pediatric Kit contains: 1 each P1753 zRIP driver module, 2 each P1133 zRIP wiresets, 2 each 1075970 zRIP DuraBelts, pediatric	1076337

Description	Part number
zRIP DuraBelt startup kits (continued)	
Other (continued)	
zRIP DuraBelt kit, adult, QDC Kit contains: 1 each P1176 abdomen cable, 1 each P1177 ribcage cable and 2 each 1073622 zRIP DuraBelts, adult	1075968
zRIP DuraBelt kit, pediatric, QDC Kit contains: 1 each P1176 abdomen cable, 1 each P1177 ribcage cable and 2 each 1075970 zRIP DuraBelts, pediatric	1075967
zRIP modules	
zRIP module, 1.5mm connectors	P1732
zRIP SUM RT module, 1.5mm connectors	P1753
zRIP module, Alice 5, 1.5mm connectors	P1736R
zRIP DuraBelt effort sensors	
zRIP DuraBelt adult, size: 24" - 74", 1.5mm panel mount black	1073622
zRIP DuraBelt adult, large, size: 45" - 123", 1.5mm panel mount black	1075973
zRIP DuraBelt, pediatric, size: 16" - 42", 1.5mm panel mount black	1075970
zRIP wiresets and cables	
zRIP wireset, 249 connector (2-pin keyhole)	P1133
zRIP wireset, Medcare, 249 connector (2-pin keyhole)	P1137
zRIP wireset, Alice PDx, 249 connector (2-pin keyhole)	P1170
zRIP cable abdomen, proprietary QDC connector	P1176
zRIP cable ribcage, proprietary QDC connector	P1177



Respiratory inductance plethysmography effort sensors



ezRIP combination module/wireset sensor

- Low initial RIP investment
- Minimizes bedside module congestion
- Easy configurability – the driver module clips directly onto the zRIP DuraBelt and legacy zRIP belt
- Non-replaceable batteries
- One-year warranty

Description	Part number
ezRIP DuraBelt kits	
ezRIP adult kit, 1.5mm, 84"	1099514
1 Abd, 1 Thx, 2 zRIP DuraBelts, adult	
ezRIP adult kit, Alice 5/LE, 84"	1099515
1 Abd, 1 Thx, 2 zRIP DuraBelts, adult	
ezRIP adult kit, 249, 84"	1099516
1 Abd, 1 Thx, 2 zRIP DuraBelts, adult	
ezRIP adult kit, 249, 16"/20"	1099517
1 Abd, 1 Thx, 2 zRIP DuraBelts, adult	
ezRIP adult, Stardust II kit	1099558
1 Thx, 1 zRIP DuraBelt, adult	
ezRIP combination module/wireset (individual component)	
KezRIP module, abdomen, 1.5mm, 84"	1094094
ezRIP module, thorax, 1.5mm, 84"	1094096
ezRIP module, abdomen, 1.5, 84", Alice 5/LE	1094110
ezRIP module, thorax, 1.5, 84", Alice 5/LE	1094111
ezRIP module, abdomen, 249, 84"	1094093
ezRIP module, thorax, 249, 84"	1094095
ezRIP module, abdomen, 249, 20"	1094920
ezRIP module, thorax, 249, 16"	1094921
ezRIP module, thorax, harwin, 26"	1094097

Understanding RIP signals vs. piezo signals

RIP	Piezo
Directly reflects breathing activity of a patient	Indirectly reflects breathing activity of a patient
Sensor wire wraps around the patient	Sensor (piezo) most commonly located near the belt buckle
Belt only needs to be snug on the patient for an accurate signal	Belt needs to be pulled taut for an accurate signal
Gain settings do not change during the night when the signal decreases	Gain settings usually need to be altered during the night to capture the signal
The signal decreasing reflects true patient breathing effort	The signal decreasing, or changing, can reflect "trapping" or artificial artifact
Thoracic belt placement is designed to be more precise than piezo	Belt placement can be anywhere in the thoracic area and pulled tight

Piezo respiratory effort sensors

Description	Part number
Piezo respiratory effort kits	
Piezo 1 channel adult kit	P1461R
Kit contains: 1 each P1420 adult Velcro® strap kit 1.5", 1 each P1460 piezo sensor 1.5", 1.5mm connectors	
Piezo 2 channel plus adult kit	P1463R
Kit contains: 2 each P1420 adult Velcro strap kit 1.5", 3 each P1460 piezo sensor 1.5", 1.5mm connectors	
Piezo 1 channel pediatric kit	P1483
Kit contains: 1 each P1416 pediatric Velcro strap kit 1.5", 1 each P1460 piezo sensor 1.5", 1.5mm connectors	
Piezo 1 channel pediatric kit	P1462R
Kit contains: 1 each P1448 pediatric Velcro strap kit 1", 1 each P1453 pediatric piezo sensor 1", 1.5mm connectors	
Piezo 2 channel plus pediatric kit	P1464R
Kit contains: 2 each P1448 pediatric Velcro strap kit 1", 3 each P1453 pediatric piezo sensor 1", 1.5mm connectors	
Adult Velcro strap kit, 1.5"	P1420
Kit contains: 1 each adult Velcro strap 1.5", size 1, size 2, and size 3	
Pediatric Velcro strap kit, 1.5"	P1416
Kit contains: 1 each pediatric Velcro strap 1.5", size 1, size 2, and size 3	
Pediatric Velcro strap kit 1"	P1448
Kit contains: 1 each pediatric Velcro strap 1", size 1, size 2, size 3, and size 4	
Piezo respiratory effort sensors	
Piezo sensor, 1.5", adult, 1.5mm connectors	P1460
Piezo sensor, 1" pediatric, 1.5mm connectors	P1453
CT2 piezo sensor, adult, 1.5mm connectors	P1582
CT2 piezo sensor, pediatric, 1.5mm connectors	P1585
Stardust CT2 effort sensor, adult, harwin connector	P1805
CT2 effort sensor, adult, Alice 5, 1.5mm connectors	P1586
CT2, effort sensor, pediatric, Alice 5, 1.5mm connectors	P1587
CT2, effort sensor, adult, Alice 4, LEMO triaxial connector	P1804
CT2, effort sensor, adult, Embletta, 249 connector (2-pin keyhole)	P1800
CT2, effort sensor, chest, adult, dark gray, 16", 249 connector (2-pin keyhole)	P1823



Piezo respiratory effort sensors

- Sensor straps are easy to size on most patients
- Low profile of the sensor straps enhances patient comfort
- Compatible with most PSG systems
- Plugs directly into standard AC channels
- Latex free
- One-year warranty



Description	Part number
Piezo respiratory effort accessories	
CT2, 18" extension strap	P1592
CT1 belt strap, adult, size 1, 1.5"w, 33.5"l	P1421
CT1 belt strap, adult, size 2, 1.5"w, 49"l	P1422
CT1 belt strap, adult, size 3, 1.5"w, 77"l	P1423

Body position sensors



Body position sensors

- Designed with patient comfort in mind
- Ability to identify the primary sleep positions (prone, supine, left, right, and upright)
- Patient's body position can be continuously monitored regardless of the amount of movement
- Attaches to the thoracic respiratory effort belt easily and comfortably
- Compatible with most PSG systems
- One-year warranty

Description	Part number
Body position sensors	
SPI (sleep position indicator) resistive, Cadwell, 8-8 modular plug	P1652
Body position sensor, DC with 3.5mm connector	P1664
Body position sensor, DC with 1.5mm connector	P1664-WS2074
Body position sensor, Alice 6, 249 connector (2-pin keyhole)	1078758
Body position sensor, Alice 5, 1.5mm connectors	P1694
Body position sensor, Alice 4, LEMO triaxial connector	P1707
Body position sensor, Alice 3, LEMO connector	P1706

Snore sensors and microphones



Snore sensors and microphones

- The snore sensor, with its thin piezo crystal, is designed to provide an accurate, reproducible signal
- The snore microphone is designed to provide a thin baseline and clear, responsive signals
- Compatible with most PSG systems
- Six-month warranty for snore sensors
- Ninety-day warranty for snore microphone

Description	Part number
Snore sensors and microphones	
Piezo snore sensor, 1.5mm connector	P1696
Snore sensor, Embletta, 249 connector (2-pin keyhole)	P1718
Snore microphone, 1.5mm connector	P1716
Snore microphone, Alice 6, 249 connector (2-pin keyhole)	1078753

Periodic limb movement

Periodic limb movement sensors

- Detect a patient's leg or arm movement
- Designed to provide clean, reliable signals while providing exceptional patient comfort
- Sensitive to subtle movements
- Do not require skin preparation prior to application
- Using kits, multiple sensors can be configured to combine PLM signals onto a single channel
- Six-month warranty

Description	Part number
Periodic limb movement (PLM)	
PLM sensor, Alice 6, 249 connector (2-pin keyhole)	1078752
Piezo PLM, 1.5mm connectors	P1697
Piezo PLM kit	P1698R
Kit contains: 1 each P1697 PLM sensor, 1 each P1689 Velcro strap PLM large, 1 each P1690 Velcro strap PLM small	
Piezo PLM 2 per kit	P1699R
Kit contains: 2 each P1697 PLM sensor, 2 each P1689 Velcro strap PLM large, 2 each P1690 Velcro strap PLM small, 1 each, P1097 safety male jumper wireset	

Periodic limb movement (PLM) accessories

Jumper, wireset, 1.5mm	P1097
Velcro strap PLM, large	P1689
Velcro strap PLM, small	P1690



Appendix: AASM scoring manual guidelines and Pro-Tech sensor application

AASM scoring manual (Version 2.3) guideline – adults	Pro-Tech sensor
Thermal airflow Adults VIII.1.A: “The sensor to detect absence of airflow...is an oronasal thermal sensor.”	Oronasal thermocouple P1222
Flow nasal pressure Adults VIII.1.B: “The sensor for detection of... hypopnea is a nasal air pressure transducer...”	PTAF2 startup kit P1287 Nasal cannulas: adult (60 per case) P1259
Respiratory effort Adults VIII.1.C: “The sensor for detection of respiratory effort is.... inductance plethysmography.”	Respiratory Inductance Plethysmography: zRIP DuraBelt kit, adult 1076329 or ezRIP adult kit, 1099514
AASM scoring manual (Version 2.3) guideline – children	Pro-Tech sensor
Thermal airflow Children VIII.1.A: “The sensor to detect absence of airflow...is an oronasal thermal sensor.”	Oronasal thermocouple P1223
Flow nasal pressure Children VIII.1.B: “The sensor for detection of... hypopnea is a nasal air pressure transducer...”	PTAF2 startup kit P1287 Nasal cannulas: pediatric (60 per case) P1263
Respiratory effort Children VIII.1.C: “The sensor for detection of respiratory effort is.... inductance plethysmography.”	Respiratory Inductance Plethysmography: zRIP DuraBelt kit, pediatric 1076330 or ezRIP module, abdomen, 1094094, ezRIP module, thorax, 1094096 and zRIP DuraBelt, pediatric, 1075970
Hypoventilation Children VIII.1.E: “Acceptable methods for assessing alveolar hypoventilation.... end-tidal PCO ₂ monitoring.” Enable EtCO ₂ and pressure monitoring simultaneously from both nostrils using one cannula.	Pro-Flow MultiPurpose Cannula: pediatric (30 per case) P1346





Medical education

For respiratory therapists, sleep technicians and nurses, Philips Respironics CEU-accredited programs target noninvasive and invasive ventilation, oxygen, and sleep therapy, and are focused on your continuing education requirements.

To learn about additional training opportunities relating to sleep apnea, visit our Partners in Training web site at <http://www.ganESCO.com/philips-attendeE/login.php>.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.

CEU programs

Sleep

The Mask Experience

1 CEU, AARC, AAST, CA BRN

A lecture to learn how to select the best mask to improve patient experience and adherence to the therapy.

- Describe the challenges the OSA patient has to be compliant to his or her CPAP therapy
- Describe why your patient needs to be compliant to CPAP therapy
- Identify the types of masks, and understand their benefits and drawbacks
- Identify how to choose the right mask to help your patient be compliant to CPAP therapy
- Discuss the significance of resupply
- Discuss recent research evidence

Fundamentals of Obstructive Sleep Apnea

1 CEU, AARC, AAST, CCMC, CA BRN

Understand the pathology and prevalence of OSA and learn how to recommend appropriate diagnostic and treatment options

- Restate the types of Sleep Disordered Breathing (SDB)
- Recognize the prevalence of SDB
- Identify and recommend appropriate diagnosis and screening procedures
- Discuss various treatment and follow-up management options.

Advanced Algorithms for SDB

1.5 CEU, AARC, AAST

Learn about leak and triggering technology, understand the CPAP breathing algorithms and identify the breathing algorithms for periodic breathing and respiratory insufficiency

- Restate the technology of auto adjusting expiratory pressure, servo ventilation, volume assured pressure support ventilation, automatic adjustment features, and comfort features
- Describe the disease states that the algorithm should be used for

Titration Guidelines for CPAP, Bi-level, Servo Ventilation and volume assured pressure support ventilation

1.5 CEU, AARC, AAST

List the titration goals for OSA management, understand the suggested protocols for titrations and review the titration process.

- Describe a titration procedure for basic CPAP
- Describe a titration procedure for servo ventilation
- Describe a titration procedure for volume assured pressure support ventilation.

Respiratory

Maximizing PAP Therapy through Advanced Techniques: servoventilation and volume assured pressure support ventilation

1 CEU, AARC, AAST

This course provides information to be able to identify the proper disease state for the use of servo ventilation or volume assured pressure support.

- Identify the link between OSA and CVD or OHS
- Describe the prevalence of OSA and CVD or OHS
- Describe the advanced techniques for treatment of OSA and CVD or OHS patients

Servo ventilation for SDB

1 CEU, AARC, AAST

This will cover the current disease processes that can be treated with servoventilation

- Recognize the breathing patterns treated with servoventilation
- Identify the disease process that the breathing patterns are seen in.
- Discuss the clinical studies that back up the servoventilation treatment for these diseases.

Home Ventilation Workshop

2 CEU, AARC, CA BRN

Learn more about the operation of the Trilogy100 and Trilogy200 life support ventilators.

- Describe the circuit options, breath types, modes, alarms, features and troubleshooting.
- Review the unique mouthpiece ventilation mode.

Mechanical insufflation-exsufflation (MI-E) technique Workshop

1 CEU, AARC

This workshop educates respiratory therapists on mechanical insufflation-exsufflation (MI-E). The session will cover indication, contraindications, operation procedure recommended treatment practice and review of a clinical case study.

- Describe the indications, contraindications and warnings/cautions for the use of mechanical insufflation-exsufflation (MI-E) techniques.
- Identify the advantages of the noninvasive use of mechanical insufflation-exsufflation (MI-E) technique over invasive airway clearance.
- Restate the details of a clinical case study
- Discuss the recommended treatment plan.





CEU programs (continued)

Respiratory (continued)

COPD Management workshop: Reducing Readmissions

- Define Triple Aim Ideology in healthcare.
- Define the term “readmission”.
- Discuss conditions typically associated with increased readmission rates.
- Describe the key factors in COPD readmission.
- Discuss the IDEAL discharge system and General Assessment of Preparedness.
- Describe intervention and implementation of a reduced readmission program using examples.
- Explain Plan-Do-Check-Act process with related to reduced readmission.
- Review and develop end-to-end treatment solutions (according to physician order) for the COPD patient “new to the system” and “chronic and considered recidivism for COPD hospitalization”.

COPD Management workshop: Part 2

- Discuss communication and education of a patient
- Recognize the common barriers found in adult learning and discuss strategies to overcome them
- Discuss resources available for COPD Patients
- Discuss the different kind of medicines commonly used in COPD management
- Discuss the required elements of personalized medication care for COPD patients*

On line only

Bi-level Rescue Program

1 CEU, AARC, AAST

Self-directed monograph highlights bi-level therapy in OSA, bi-level rescue study results, implementation and protocols.

- Describe the use of bi-level therapy as an option for patients

An Introduction to Sleep Diagnostic Equipment

3 CEU's, AAST

Self-directed monograph with 3 chapters on nasal pressure flow systems, pulse transit time, and respiratory inductance plethysmography

- Identify the difference between various nasal pressure flow systems, pulse transit time, and respiratory inductance plethysmography

Mask Monograph

1 CEU, AARC, AAST

Self-directed monograph that focuses on choosing the right mask for the right patient

- Discuss how to choose the right mask to help your patient be compliant to CPAP therapy
- Describe the different types of masks available and when they should be used

Motivational Enhancement therapy (MET)

1 CEU, AARC, AAST

Review MET as it relates to improving adherence to PAP treatment

- Describe the specific behavioral interventions and theories about behavioral change
- Recognize when to intervene.

On Line Monographs at www.ganESCO.com



Value-added programs

Philips Respironics offers a broad array of programs aimed at helping to maximize your growth and operational efficiencies. From compliance and resupply to training and warranties, we provide a spectrum of unique programs focused on facilitating success for you and your patients.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.

Device programs

Flex Promise program

Our Flex Promise Program is Philips Respironics' commitment that if any patient is noncompliant on any of our Flex technologies, upon return of the unit with Flex technology, we will apply a 100% credit of the corresponding unit purchase price toward the purchase of a Philips Respironics bi-level device with Bi-Flex.

Bi-level Rescue program

This program is intended to help patients continue their PAP therapy when they are not able to tolerate standard treatment.

Extended Warranty program

The Extended Warranty supports all DreamStation devices and adds three years of coverage to the standard two-year warranty – providing five years of protection on a DreamStation device or humidifier against unexpected repair costs and shortening the time the device is out of service.

DreamStation Non-warranty Service program

For devices not covered under warranty, the DreamStation Non-warranty Service Program offers a simple two-tiered approach to repair costs that helps to accelerate the return of the repaired unit to the patient and reimbursement turnaround times.

Mask programs

30-Day Mask Satisfaction Promise program

Your patient will be satisfied with our mask or we'll replace it. If a patient discontinues use of a program-approved mask in favor of another mask during the first 30 days of use, we will replace the original mask at no charge. Certain terms and conditions apply.

Ventilation programs

Ventilation Promise program

Progressive diseases like COPD and ALS have one thing in common – a patient's ventilator needs will change over time. We want to make it easy for you to transition your patients to higher level ventilator devices as their needs become more challenging. To enable you to easily transition patients to other devices as their diseases progress, the Ventilation Promise Program will provide a credit equal to the purchase price of any Philips Respironics' System One device to be applied toward the purchase of a BiPAP AVAPS or Trilogy ventilator. Certain terms and conditions apply.

One COPD program

One COPD takes a patient-centric approach—based on clinical research and backed by Philips devices, support, and expertise. Technology, education, and connectivity all work together to provide a seamless care experience geared towards patient satisfaction, better patient care, and reduced cost of care.

Customer programs

Cognitive Behavioral Therapy for Insomnia (CBTI)

This program provides sleep labs with a means to implement a behavioral sleep medicine program or enhance their current offerings without incurring the development cost and time of doing it themselves. Additionally, it enables them to implement an objective tool (Actiwatch) for therapeutic assessment as an adjunct to traditional, subjective diaries.

Know How Webinars

The Philips Respironics Know How program is a series of live and on-demand webinars that help to address key business and clinical challenges and opportunities faced by homecare providers and healthcare professionals. Visit <http://www.learningconnection.philips.com/en/sleep-and-respiratory-care-education> to learn more and to register for the various sessions.



