


Philips NatureDynamics Mounting Clip

General Information


Universal mounting clip for NatureDynamics tubular lighting. The clip has provisions for securing the orientation of the tube.

Mounting options

Mounted on 1 steel wire



Mounting on 2 steel wires



Mounted by screw



Orientation lock feature



Product Data

Order product name:
Philips NatureDynamics Mounting Clip
Numerator - Quantity Per Pack 100
Material Number: (12NC) 912500100197
Net Weight (Piece) 0.016 kg

Dimensions

L x W x H: 41.5 x 50.1 x 36.7 mm
Suspension wire max diameter: 4mm
Max screw diameter: 4mm
Center distance for 2 suspension wire option: 32mm
Material: PA