



1.0 GS1 Barcode Format



- 1 Application Identifier, AI (01), indicates that the GS1 data field contains a GTIN
- 2 GTIN-14
- 3 AI (11) indicates that the GS1 AI data field contains date of manufacture (applicable to MDR compliant labels)
- 4 Date of manufacture (YYMMDD, example: 210914) (applicable to MDR compliant labels)
- 5 AI (17) indicates that the GS1 AI data field contains Expiration Date
- 6 Expiration Date (YYMMDD, example: 230914)
- 7 AI (10) indicates that the GS1 AI data field contains Lot Number
- 8 Lot Number (XXXXXXXXXX, 9 characters)

2.0 GTIN-14 Barcode Format


















- 1 indicator digit for different levels of product hierarchy (for product with package of 1 unit per box: the 0 indicator is for both pouch and box; for product with package of multiple units: 0 is for pouch and 2 is for box)
- 2 GS1 Company Prefix (SPNC Prefix: 081313202)
- 3 Item Reference Number (example: 074)
- 4 Calculated check digit reference (calculation based on previous 13 digits to ensure global uniqueness)

3.0 AngioSculpt PTCA GTIN Number and Barcode

Model Number	Configuration	GTIN Number	GTIN Portion of the GS1 Barcode
2027-2010	2.0 mm x 10 mm, OTW, PTCA	00813132021016	 00813132021016
2027-2015	2.0 mm x 15mm, OTW, PTCA	00813132021023	 00813132021023
2027-2020	2.0 mm x 20 mm, OTW, PTCA	00813132021030	 00813132021030
2027-2510	2.5 mm x 10 mm, OTW, PTCA	00813132021047	 00813132021047
2027-2515	2.5 mm x 15 mm, OTW, PTCA	00813132021054	 00813132021054
2027-2520	2.5 mm x 20 mm, OTW, PTCA	00813132021061	 00813132021061
2027-3010	3.0 mm x 10 mm, OTW, PTCA	00813132021078	 00813132021078
2027-3015	3.0 mm x 15 mm, OTW, PTCA	00813132021085	 00813132021085
2027-3020	3.0 mm x 20 mm, OTW, PTCA	00813132021092	 00813132021092
2027-3510	3.5 mm x 10 mm, OTW, PTCA	00813132021108	 00813132021108







Model Number	Configuration	GTIN Number	GTIN Portion of the GS1 Barcode
2027-3515	3.5 mm x 15 mm, OTW, PTCA	00813132021115	 00813132021115
2027-3520	3.5 mm x 20 mm, OTW, PTCA	00813132021122	 00813132021122
2200-2006	2.0 mm x 6 mm, RX, PTCA	00813132021139	 00813132021139
2200-2010	2.0 mm x 10 mm, RX, PTCA	00813132021146	 00813132021146
2200-2015	2.0 mm x 15 mm, RX, PTCA	00813132021153	 00813132021153
2200-2506	2.5 mm x 6 mm, RX, PTCA	00813132021160	 00813132021160
2200-2510	2.5 mm x 10 mm, RX, PTCA	00813132021177	 00813132021177
2200-2515	2.5 mm x 15 mm, RX, PTCA	00813132021184	 00813132021184
2200-3006	3.0 mm x 6 mm, RX, PTCA	00813132021191	 00813132021191
2200-3010	3.0 mm x 10 mm, RX, PTCA	00813132021207	 00813132021207

Model Number	Configuration	GTIN Number	GTIN Portion of the GS1 Barcode
2200-3015	3.0 mm x 15 mm, RX, PTCA	00813132021214	 00813132021214
2200-3506	3.5 mm x 6 mm, RX, PTCA	00813132021221	 00813132021221
2200-3510	3.5 mm x 10 mm, RX, PTCA	00813132021238	 00813132021238
2200-3515	3.5 mm x 15 mm, RX, PTCA	00813132021245	 00813132021245
2001-2010	2.0 mm x 10 mm, RX, PTCA	00813132022075	 00813132022075
2001-2015	2.0 mm x 15 mm, RX, PTCA	00813132022082	 00813132022082
2001-2020	2.0 mm x 20 mm, RX, PTCA	00813132022099	 00813132022099
2001-2510	2.5 mm x 10 mm, RX, PTCA	00813132022105	 00813132022105
2001-2515	2.5 mm x 15 mm, RX, PTCA	00813132022112	 00813132022112
2001-2520	2.5 mm x 20 mm, RX, PTCA	00813132022129	 00813132022129
2001-3010	3.0 mm x 10 mm, RX, PTCA	00813132022136	 00813132022136

Model Number	Configuration	GTIN Number	GTIN Portion of the GS1 Barcode
2001-3015	3.0 mm x 15 mm, RX, PTCA	00813132022143	 00813132022143
2001-3020	3.0 mm x 20 mm, RX, PTCA	00813132022150	 00813132022150
2001-3510	3.5 mm x 10 mm, RX, PTCA	00813132022167	 00813132022167
2001-3515	3.5 mm x 15 mm, RX, PTCA	00813132022174	 00813132022174
2001-3520	3.5 mm x 20 mm, RX, PTCA	00813132022181	 00813132022181

4.0 AngioSculpt Evo GTIN Number and Barcode

Model Number	Configuration	GTIN Number	GTIN Portion of the GS1 Barcode
2200-2006-B	2.0 mm x 06 mm, Evo RX PTCA	00813132026790	 00813132026790
2200-2010-B	2.0 mm x 10 mm, Evo RX PTCA	00813132026806	 00813132026806
2200-2015-B	2.0 mm x 15 mm, Evo RX PTCA	00813132026813	 00813132026813
2200-2020-B	2.0 mm x 20 mm, Evo RX PTCA	00819192026820	 00819192026820
2200-2506-B	2.5 mm x 06 mm, Evo RX PTCA	00813132026837	 00813132026837
2200-2510-B	2.5 mm x 10 mm, Evo RX PTCA	00813132026844	 00813132026844
2200-2515-B	2.5 mm x 15 mm, Evo RX PTCA	00813132026851	 00813132026851
2200-2520-B	2.5 mm x 20 mm, Evo RX PTCA	00813132026868	 00813132026868
2200-3006-B	3.0 mm x 06 mm, Evo RX PTCA	00813132026875	 00813132026875

Model Number	Configuration	GTIN Number	GTIN Portion of the GS1 Barcode
2200-3010-B	3.0 mm x 10 mm, Evo RX PTCA	00813132026882	 00813132026882
2200-3015-B	3.0 mm x 15 mm, Evo RX PTCA	00813132026899	 00813132026899
2200-3020-B	3.0 mm x 20 mm, Evo RX PTCA	00813132026905	 00813132026905
2200-3506-B	3.5 mm x 06 mm, Evo RX PTCA	00813132026912	 00813132026912
2200-3510-B	3.5 mm x 10 mm, Evo RX PTCA	00813132026929	 00813132026929
2200-3515-B	3.5 mm x 15 mm, Evo RX PTCA	00813132026936	 00813132026936
2200-3520-B	3.5 mm x 20 mm, Evo RX PTCA	00813132026943	 00813132026943