


## GS1 / HIBC BAR CODE SHEET

### GS1-128 Barcode Format



- 1 Application Identifier, AI (01), indicates that the GS1 AI data field contains a GTIN
- 2 GTIN-14
- 3 AI (17) indicates that the GS1 AI data field contains expiration date
- 4 Expiration Date (YYMMDD, example: 160910)
- 5 AI(10) indicates that the GS1 AI data field contains Lot number
- 6 Lot number (XXXXXXXXXX, 9 characters)

### GTIN-14 Barcode Format







- 1 Indicator Digit for different levels of product hierarchy (For product with package of 1 unit per box: the 0 indicator is for both pouch and box; For product with package of multiple units: 0 is for pouch and 2 is for box)
- 2 GS1 Company Prefix (SPNC Prefix: 081313202)
- 3 Item Reference Number (example: 074)
- 4 Calculated Check Digit Reference (calculation based on previous 13 digits to ensure global uniqueness).



# QuickCat™

Extraction Catheter

## GS1 / HIBC BAR CODE SHEET

| Model #  | Configuration | GTIN Portion of the GS1 Barcode   | HIBC Barcode   |
|----------|---------------|---|--|
|          |               | <b>GTIN-14 BOX BARCODE (Updated)</b>  | <b>HIBC BOX BARCODE (Original)</b>   |
| 60090-01 | 6F            | <br>00813132020989 | <br>+M2046009001AA1K  |
|          |               | <b>GTIN-14 POUCH BARCODE (Updated)</b>  | <b>HIBC POUCH BARCODE (Original)</b>   |
| 60090-01 | 6F            | <br>00813132020989 | <br>+ M2046009001AA0J |