

Turbo-Elite™


Laser Atherectomy Catheter GS1 BARCODE SHEET

GS1-128 Barcode Format



- 1 Application Identifier, AI (01), indicates that the GS1 AI data field contains a GTIN
- 2 GTIN-14
- 3 AI (17) indicates that the GS1 AI data field contains expiration date
- 4 Expiration Date (YYMMDD, example: 160910)
- 5 AI(10) indicates that the GS1 AI data field contains Lot number
- 6 Lot number (XXXXXXXXXX, 9 characters)


















GTIN-14 Barcode Format



- 1 Indicator Digit for different levels of product hierarchy (For product with package of 1 unit per box: the 0 indicator is for both pouch and box; For product with package of multiple units: 0 is for pouch and 2 is for box)
- 2 GS1 Company Prefix (SPNC Prefix: 081313202)
- 3 Item Reference Number (example: 074)
- 4 Calculated Check Digit Reference (calculation based on previous 13 digits to ensure global uniqueness).

Turbo-Elite™

Laser Atherectomy Catheter GS1 BARCODE SHEET

Model #	Configuration	GTIN Portion of the GS1 Barcode	GTIN Portion of the GS1 Barcode
		GTIN-14 BOX BARCODE (Updated)	GTIN-14 BOX BARCODE (Previous)
425-135-02 (OUS)	2.5 mm catheter TER OTW, 0.035 in. compatible	 00813132024826	 00813132020224
425-135-01 (US)	2.5 mm catheter TER OTW, 0.035 in. compatible	 00813132024819	 00813132020217
423-135-02 (OUS)	2.3 mm catheter TER OTW, 0.035 in. compatible	 00813132024802	 00813132020200
423-135-01 (US)	2.3 mm catheter TER OTW, 0.035 in. compatible	 00813132024796	 00813132020194
425-011	2.5 mm catheter OTW 0.018 in. compatible	 00813132024789	 00813132020187
423-001	2.3 mm catheter OTW 0.018 in. compatible	 00813132024772	 00813132020170
420-159	2.0 mm catheter RX 0.014 in. compatible	 00813132024765	 00813132020163
420-006	2.0 mm catheter OTW 0.018 in. compatible	 00813132024758	 00813132020156
417-156	1.7 mm catheter RX 0.014 in. compatible	 00813132024741	 00813132020149

The Spectranetics Corporation











9965 Federal Drive, Colorado Springs, CO 80921

Tel: 719-447-2000 • Fax: 719-447-2022

Implemented: 07Mar2017














Turbo-Elite™

Laser Atherectomy Catheter GS1 BARCODE SHEET

Model #	Configuration	GTIN Portion of the GS1 Barcode	GTIN Portion of the GS1 Barcode
		GTIN-14 BOX BARCODE (Updated)	GTIN-14 BOX BARCODE (Previous)
417-152	1.7 mm catheter OTW 0.018 in. compatible	 00813132024734	 00813132020132
414-159	1.4 mm catheter RX 0.014 in. compatible	 00813132024727	 00813132020125
414-151	1.4 mm catheter OTW 0.014 in. compatible	 00813132024710	 00813132020118
410-154	0.9 mm catheter RX 0.014 in. compatible	 00813132024703	 00813132020101
410-152	0.9 mm catheter OTW 0.014 in. compatible	 00813132024697	 00813132020095











Turbo-Elite™

Laser Atherectomy Catheter GS1 BARCODE SHEET

Model #	Configuration	GTIN Portion of the GS1 Barcode	GTIN Portion of the GS1 Barcode
		GTIN-14 POUCH BARCODE (Updated)	GTIN-14 POUCH BARCODE (Previous)
425-135-02 (OUS)	2.5 mm catheter TER OTW, 0.035 in. compatible	 00813132024826	 00813132020224
425-135-01 (US)	2.5 mm catheter TER OTW, 0.035 in. compatible	 00813132024819	 00813132020217
423-135-02 (OUS)	2.3 mm catheter TER OTW, 0.035 in. compatible	 00813132024802	 00813132020200
423-135-01 (US)	2.3 mm catheter TER OTW, 0.035 in. compatible	 00813132024796	 00813132020194
425-011	2.5 mm catheter OTW 0.018 in. compatible	 00813132024789	 00813132020187
423-001	2.3 mm catheter OTW 0.018 in. compatible	 00813132024772	 00813132020170
420-159	2.0 mm catheter RX 0.014 in. compatible	 00813132024765	 00813132020163
420-006	2.0 mm catheter OTW 0.018 in. compatible	 00813132024758	 00813132020156
417-156	1.7 mm catheter RX 0.014 in. compatible	 00813132024741	 00813132020149

Turbo-Elite™

Laser Atherectomy Catheter GS1 BARCODE SHEET

Model #	Configuration	GTIN Portion of the GS1 Barcode	GTIN Portion of the GS1 Barcode
		GTIN-14 POUCH BARCODE (Updated)	GTIN-14 POUCH BARCODE (Previous)
417-152	1.7 mm catheter OTW 0.018 in. compatible	 00813132024734	 00813132020132
414-159	1.4 mm catheter RX 0.014 in. compatible	 00813132024727	 00813132020125
414-151	1.4 mm catheter OTW 0.014 in. compatible	 00813132024710	 00813132020118
410-154	0.9 mm catheter RX 0.014 in. compatible	 00813132024703	 00813132020101
410-152	0.9 mm catheter OTW 0.014 in. compatible	 00813132024697	 00813132020095