


Quick-Access™


Needle Holder

GS1 AND HIBC BARCODE SHEET

GS1-128 Barcode Format




(01)00813132020743(17)160910(10)XXXXXXXXXX




- 1 Application Identifier, AI (01), indicates that the GS1 AI data field contains a GTIN
- 2 GTIN-14
- 3 AI (17) indicates that the GS1 AI data field contains expiration date
- 4 Expiration Date (YYMMDD, example: 160910)
- 5 AI(10) indicates that the GS1 AI data field contains Lot number
- 6 Lot number (XXXXXXXXXX, 9 characters)

GTIN-14 Barcode Format











00813132020743



- 1 Indicator Digit for different levels of product hierarchy (For product with package of 1 unit per box: the 0 indicator is for both pouch and box; For product with package of multiple units: 0 is for pouch and 2 is for box)
- 2 GS1 Company Prefix (SPNC Prefix: 081313202)
- 3 Item Reference Number (example: 074)
- 4 Calculated Check Digit Reference (calculation based on previous 13 digits to ensure global uniqueness).

Quick-Access™
Needle Holder

GS1 AND HIBC BARCODE SHEET

Model #	Configuration	GTIN Portion of the GS1 Barcode	HIBC Barcode
		GTIN-14 BOX BARCODE (Updated)	HIBC BOX BARCODE (Original)
519-001	Quick-Access Needle Holder (Single Pack)	 00813132020729	 +M204519001AA1K
519-005	Quick-Access Needle Holder (Five Pack)	 20813132020730	 +M204519005AA1O
		GTIN-14 POUCH BARCODE (Updated)	HIBC POUCH BARCODE (Original)
519-001	Quick-Access Needle Holder (Single Pack)	 00813132020729	 +M204519001AA0J
519-005	Quick-Access Needle Holder (Five Pack)	 00813132020736	 +M204519005AA0N