



Turbo-Booster™

Laser Guide Catheter



GS1 AND HIBC BARCODE SHEET

GS1-128 Barcode Format



- 1 Application Identifier, AI (01), indicates that the GS1 AI data field contains a GTIN
- 2 GTIN-14
- 3 AI (17) indicates that the GS1 AI data field contains expiration date
- 4 Expiration Date (YYMMDD, example: 160910)
- 5 AI(10) indicates that the GS1 AI data field contains Lot number
- 6 Lot number (XXXXXXXXXX, 9 characters)

GTIN-14 Barcode Format











- 1 Indicator Digit for different levels of product hierarchy (For product with package of 1 unit per box: the 0 indicator is for both pouch and box; For product with package of multiple units: 0 is for pouch and 2 is for box)
- 2 GS1 Company Prefix (SPNC Prefix: 081313202)
- 3 Item Reference Number (example: 074)
- 4 Calculated Check Digit Reference (calculation based on previous 13 digits to ensure global uniqueness).

Turbo-Booster™

Laser Guide Catheter

GS1 AND HIBC BARCODE SHEET

Model #	Configuration	GTIN Portion of the GS1 Barcode	HIBC Barcode
		GTIN-14 BOX BARCODE (Updated)	HIBC BOX BARCODE (Original)
518-063	8 Fr Turbo-Booster 3-Pk	 20813132020280	 +M204518063AA1R
518-043	7 Fr Turbo-Booster 3-Pk	 20813132020273	 +M204518043AA1P
		GTIN-14 POUCH BARCODE (Updated)	HIBC POUCH BARCODE (Original)
518-063	8 Fr Turbo-Booster 3-Pk	 00813132020286	 +M204518063AA0Q
518-043	7 Fr Turbo-Booster 3-Pk	 00813132020279	 +M204518043AA0O