





Spectranetics®

Turbo-Booster®

Laser Guide Catheter

BAR CODE SHEET

Model #	Configuration	Barcode #	Bar Code
BOX BARCODES			
518-043	7F	+M20451804326	 +M20451804326\
518-063	8F	+M20451806328	 +M20451806328#
POUCH BARCODES			
518-043	7F	+M20451804315	 +M20451804315C
518-063	8F	+M20451806317	 +M20451806317q

© 2009 Spectranetics Approved for External Distribution D007817-00 052009

9965 Federal Dr., Colorado Springs, Colorado 80921 • 719-447-2000 • Fax: 719-447-2022 • Customer Service: 800-231-0978 • www.spnc.com

implemented: 13May09