











Turbo-Elite™

Laser Atherectomy Catheter

BAR CODE SHEET




Model #	Configuration	Barcode #	Bar Code
BOX BARCODES			
425-135-02 (OUS)	2.5 mm catheter TER OTW, 0.35 in. compatible	+M2044231350216	 +M2044231350216
425-135-01 (US)	2.5 mm catheter TER OTW, 0.35 in. compatible	+M2044251350115	 +M2044251350115
423-135-02 (OUS)	2.3 mm catheter TER OTW, 0.35 in. compatible	+M2044231350214	 +M2044231350214
423-135-01 (US)	2.3 mm catheter TER OTW, 0.35 in. compatible	+M2044231350113	 +M2044231350113
425-011	2.5 mm catheter OTW	+M2044250111/	 +M2044250111/
423-001	2.3 mm catheter OTW	+M2044230011.	 +M2044230011.
420-159	2.0 mm catheter RX	+M20442015915	 +M20442015915
420-006	2.0 mm catheter OTW	+M2044200061\$	 +M2044200061\$
417-156	1.7 mm catheter RX	+M20441715618	 +M20441715618
417-152	1.7 mm catheter OTW	+M20441715214	 +M20441715214
414-159	1.4 mm catheter RX	+M20441415918	 +M20441415918



Turbo-Elite™

Laser Atherectomy Catheter

BAR CODE SHEET

Model #	Model #	Model #	Model #
BOX BARCODES			
414-151	1.4 mm catheter OTW	+M20441415110	 +M20441415110
410-154	0.9 mm catheter RX	+M2044101541%	 +M2044101541%
410-152	0.9 mm catheter OTW	+M2044101521/	 +M2044101521/

Turbo-Elite™

Laser Atherectomy Catheter

BAR CODE SHEET





Model #	Configuration	Barcode #	Bar Code
POUCH BARCODES			
425-135-02 (OUS)	2.5 mm catheter TER OTW, 0.35 in. compatible	+M2044251350205	 +M2044251350205
425-135-01 (US)	2.5 mm catheter TER OTW, 0.35 in. compatible	+M2044251350104	 +M2044251350104
423-135-02 (OUS)	2.3 mm catheter TER OTW, 0.35 in. compatible	+M2044231350203	 +M2044231350203
423-135-01 (US)	2.3 mm catheter TER OTW, 0.35 in. compatible	+M2044231350102	 +M2044231350102
425-011	2.5 mm catheter OTW	+M2044250110\$	 +M2044250110\$
423-001	2.3 mm catheter OTW	+M2044230010-	 +M2044230010-
420-159	2.0 mm catheter RX	+M20442015904	 +M20442015904
420-006	2.0 mm catheter OTW	+M2044200060	 +M2044200060
417-156	1.7 mm catheter RX	+M20441715607	 +M20441715607
417-152	1.7 mm catheter OTW	+M20441715203	 +M20441715203



Turbo-Elite™

Laser Atherectomy Catheter

BAR CODE SHEET

Model #	Model #	Model #	Model #
POUCH BARCODES			
414-159	1.4 mm catheter RX	+M20441415907	 +M20441415907
414-151	1.4 mm catheter OTW	+M2044141510%	 +M2044141510%
410-154	0.9 mm catheter RX	+M2044101540+	 +M2044101540+
410-152	0.9 mm catheter OTW	+M2044101520\$	 +M2044101520\$